CerTest Rapid Test

CerTest
BIOTEC

CerTest *Clostridium difficile* antigen GDH

ONE STEP CLOSTRIDIUM DIFFICILE CARD TEST

CE IVD

The test should be used by professional in IVD devices

These instructions for use apply to the following references:

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GD820001V; GD820001VC; GD820001VCC

Initial part of the Instructions for use is common for all the references/

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Specific parts for each of the references are described in the corresponding annexes (see tables below)/

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KIT 20 CARD CLOSTRIDIUM DIFFICILE ANTIGEN GDH WITH POSITIVE CONTROL (SEE ANNEX 2) /

KIT DE 20 CARCASAS DE CLOSTRIDIUM DIFFICILE ANTIGEN GDH CON CONTROL POSITIVO (VER ANEXO 2)

PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
CerTest <i>Clostridium difficile</i> antigen GDH Card (20 tests kit with stool vials + PC)	GD820001VC

Table A 2. Reference for kit 20 card of CerTest *Clostridium difficile* antigen GDH with positive control / Referencia para kit de 20 carcassas de CerTest *Clostridium difficile* antigen GDH con control positivo.

KIT 20 CARD CLOSTRIDIUM DIFFICILE ANTIGEN GDH WITH POSITIVE AND NEGATIVE CONTROL (SEE ANNEX 3) /

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ENGLISH

2. Intended use	7
3. Summary and Explanation	7
4. Test principle	8
5. Storage and stability	8
6. Precautions for users	8
7. Specimen collection and preparation	10
8. Limitations of the test	10
9. Expected values	11
10. Performance characteristics	11
A. Analytical sensitivity (detection limit)	11
B. Clinical sensitivity and specificity	12
C. Hook Effect	13
D. Cross Reactivity	14
E. Interferences	15
F. Repeatability and Reproducibility	16
ANNEX 1	17
I. Materials	17
II. Test procedures (samples)	17
III. Interpretation of results	18
IV. Quality Control	19
ANNEX 2	20
I. Materials	20
II. Test procedures (samples)	20
III. Interpretation of results	21
IV. Quality Control	22
A. External control. CerTest C. <i>difficile</i> antigen GDH Positive Control	22
A.1 Positive control procedure	22
A.1.1 Interpretation positive control procedure	23
ANNEX 3	24

I.	Materials.....	24
II.	Test procedures (samples)	24
III.	Interpretation of results	25
IV.	Quality Control	26
	A. External positive control. CerTest C. difficile antigen GDH Positive Control	26
A.1	Positive control procedure.....	26
A.1.1	Interpretation positive control procedure	27
	B. External negative control. CerTest Negative Control.....	28
B.1	Negative control procedure	28
B.1.1.	Interpretation negative control procedure	29

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ENGLISH

2. Intended use

CerTest *Clostridium difficile* antigen GDH one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH) in stool samples.

CerTest *Clostridium difficile* antigen GDH card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Clostridium difficile* infection.

3. Summary and Explanation

Clostridium difficile (*C. difficile*), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics (1).

C. difficile is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide. Although traditionally considered an infection affecting patients in health care facilities, in the last 10-20 years *C. difficile* infection is increasingly identified in the community among people with no recent hospitalizations and without prior antibiotic exposure (2–4).

Disease due to *C. difficile* develops when the bacteria is allowed to proliferate in the colon, most commonly after antibiotic use has resulted in sufficient disturbance of the intestinal microflora. Overall, the use of antibiotics is associated with a three-fold increased risk of community-associated *C. difficile* infection (CA-CDI). However, the risk of CA-CDI is largely dependent on the class of antibiotic use (1,3).

C. difficile can release two high-molecular-weight toxins, Toxin A and Toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon, and death. *Clostridium difficile* Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent marker for the microorganism (2,5).

The use of a direct faecal GDH screen, together with a faecal Toxin AB test, could improve the diagnosis of *Clostridium difficile* infection. Patients with a positive faecal GDH but negative toxin can be monitored closely while awaiting the result of a toxigenic culture (TC), that is used as the gold standard technique to determine *Clostridium difficile* infection (5).

Preventing *C. difficile* infection offers a potentially significant improvement in patient's outcomes, as well as a reduction in hospital costs and resource expenditures. Careful use of antibiotics remains a prominent public health strategy to reduce *C. difficile* infection cases (1,3).

4. Test principle

CerTest *Clostridium difficile* antigen GDH is based on the principle of a qualitative immunochromatographic assay for qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH) in human stool samples from patients suspected of *Clostridium difficile* antigen GDH infection.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against GDH and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-GDH) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigen of the diluted sample reacts with the red-coloured conjugate complex (anti-GDH monoclonal antibodies-red polystyrene microspheres) which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-GDH antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no GDH antigen presence and yet, the antigen may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-GDH antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appear. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

5. Storage and stability

The kits should be stored at 2-30°C, in these conditions can be used until the expiration date which is stated on the kit label. All the kit components are for only one use and should remain in its primary packaging until use. Do not freeze.

6. Precautions for users

- **For professional *in vitro* diagnostic use only.**
- Read the instructions for use carefully before use the kit.
- Do not use after expiration date.

- Do not use the kit if the label that seals the outer box is broken.
- Do not use product if the protective external box or the protective aluminium pouches are opened or damaged upon arrival.
- Do not use the test (and positive/negative controls if the final reference contains them) if desiccant material is not present or broken inside the aluminium pouch.
- Components provided in the kit are approved for use in **CerTest Clostridium difficile antigen GDH**. Do not mix or use with components from other product kit or lots. Do not use with other commercial kit components (from another supplier).
- All the components of the kit are for single use. Do not reuse.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves or other personal protective equipment such as goggles and mask that will be considered necessary. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous, as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples. Each sample must be correctly and unequivocally identified, in order to guarantee the correct traceability of the samples.
- A new test must be used for each sample to avoid contaminations errors.
- Clean up spills thoroughly using an appropriate disinfectant.
- Reagents contain preservatives (<0.1% sodium azide). Avoid any contact with skin or mucous membrane. In accordance with Regulation (EC) No 1907/2006 (REACH), **CerTest Clostridium difficile antigen GDH** do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration. Material Safety Data Sheet are not included with this device, nevertheless it could be consulted upon request to CerTest Biotec.
- The certificate of analysis is not included with the device, however, it could be requested to Certest Biotec in case of need.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Visual interpretation of the results by coloured lines, interpretation of the results must be done by professional user without problems in colour interpretation.
- The tests (vials, devices, positive/negative control if any) should be discarded in a proper biohazard container after testing. These containers should be discarded in accordance with local or national laws or regulations.

7. Specimen collection and preparation


Samples should be collected in a clean container and processed as soon as possible after collection. Follow proper infection control practices.

Samples can be preserved at room temperature if they are going to be analysed into the same day after taking them. The recommendation is processed the sample immediately after collection.

If this is not possible, stool samples could be stored in the refrigerator (2-8°C) for a maximum of 7 days. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the **stool collection tube** (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 50 mg) and add it to the **stool collection tube**. For liquid samples, add approx. 125 µL in the **stool collection tube** using a micropipette.
2. Close the **stool collection tube**. Shake the tube in order to assure good sample dispersion (3).

Continue immediately with the test procedure. See Annexes. Each sample vial (containing the sample) and cassette with which will be tested that sample, should be perfectly recognizable including patient's identification.

8. Limitations of the test

- 1) The test must be carried out within 2 hours after opening the sealed bag.
- 2) An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- 3) The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.

- 4) **CerTest Clostridium difficile antigen GDH** should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- 5) Positive results determine the presence of *Clostridium difficile* in faecal samples; nevertheless, it can be due to toxigenic or non-toxigenic strains of *Clostridium difficile*. A positive result should be followed up with additional laboratory techniques (toxigenic culture) to determine the strain. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- 6) A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation persist, a *Clostridium difficile* determination should be carried out, on a sample from an enrichment culture.
- 7) Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- 8) Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

9. Expected values

Clostridium difficile is associated with 95-100% cases of pseudomembranous colitis, 60-75% cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhoea cases. In addition, nearly 223.900 people in the United States required hospital care for *C. difficile* and at least 12.800 people died in 2017. Although U.S. outpatient antibiotic prescription rate decreased from 2011 to 2016, at least 30% of outpatient antibiotic prescription are estimated to be unnecessary, which highlights the need to improve outpatient prescribing (4–6).

While the development of new treatments is one of these key actions, such investments must be coupled with dedicated efforts toward preventing infections in the first place, slowing the development of resistance through better antibiotic use, and stopping the spread of resistance (6).

10. Performance characteristics

A. Analytical sensitivity (detection limit)

Detection limit value (LoD) of **CerTest Clostridium difficile antigen GDH** is 0.39 ng/mL of *Clostridium difficile* antigen GDH. There is not any external standard approved as a general and reference standard. The internal standard used for validation of this CerTest product is Glutamate dehydrogenase recombinant protein.

B. Clinical sensitivity and specificity

An external retrospective evaluation (Hospital Clínico Lozano Blesa, Zaragoza) with frozen stool samples was performed using **CerTest Clostridium difficile antigen GDH** versus another immunochromatographic (IC) test. Discordant results between the tests were analyzed by qPCR. This is the evaluation criteria. The results were as follows:

IC test: CerTest Clostridium difficile antigen GDH	Evaluation criteria		
	+	-	Total
+	86	1	87
-	3	160	163
Total	89	161	250

Table 1 CerTest Clostridium difficile antigen GDH performance versus evaluation criteria.

CerTest Clostridium difficile antigen GDH vs Evaluation criteria		
	Mean Value	95% confidence interval
Sensitivity	96.6%	90.5-99.3%
Specificity	99.4%	96.6-100.0%
PPV	98.9%	93.8-100.0%
NPV	98.2%	94.7-99.6%

Table 2 CerTest Clostridium difficile antigen GDH clinical results.

The most recent external retrospective evaluation performed to the test (with frozen stool samples) showed the following results: sensitivity 96.6% (95% confidence interval: 90.5-99.3%) and specificity: 99.4% (95% confidence interval: 96.6-100.0%).

The results showed a high sensitivity and specificity to detect *Clostridium difficile* using **CerTest Clostridium difficile antigen GDH**.

The following tables provide information about the evaluations **CerTest Clostridium difficile antigen GDH** has undergone.

Nº	Site	Sample type	Reference	Target
1	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-May 2021	Human stool samples	(17.1.1)	<i>Clostridium difficile</i> antigen GDH
2	Hospital Clínico Lozano Blesa (Zaragoza, Spain) April-June 2017	Human stool samples	(17.1.2)	<i>Clostridium difficile</i> antigen GDH
3	Louvain, (Belgium) January- December 2016	Human stool samples	(17.1.3)	<i>Clostridium difficile</i> antigen GDH
4	ECCMID April 2016, Amsterdam	Human stool samples	(17.1.4)	<i>Clostridium difficile</i> antigen GDH

Table 3 Resume of Clinical evaluations. Site, sample type, target.

N°	TP	FP	FN	TN	N	Sens. (%) CI (%)	Spec. (%) CI (%)	PPV (%) CI (%)	NPV (%) CI (%)
1	86	1	3	160	250	96.6 (90.5-99.3)	99.4 (96.6-100.0)	98.9 (93.8-100.0)	98.2 (94.7-99.6)
2	51	0	0	75	126	100.0 (93.0-100.0)	100.0 (95.2-100.0)	100.0 (93.0-100.0)	100.0 (95.2-100.0)
3	130	8	8	332	478	94.2 (88.9-97.5)	97.6 (95.4-99.0)	94.2 (88.9-97.5)	97.6 (95.4-99.0)
4	27	8	6	117	158	81.8 (64.5-93.0)	93.6 (87.8-97.2)	77.1 (59.9-89.6)	95.1 (89.7-98.2)
Total	294	17	17	684	1012	94.5 (91.4-96.8)	97.6 (96.1-98.6)	94.5 (91.4-96.8)	97.6 (96.1-98.6)

Table 4 Clinical evaluations and global results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI) for CerTest *Clostridium difficile* antigen GDH.

C. Hook Effect

CerTest *Clostridium difficile* antigen GDH does not show full inhibitory hook effect at concentration higher than $10^5 \times \text{LoD}$.

D. Cross Reactivity

An evaluation was performed to determine the cross reactivity of **CerTest Clostridium difficile antigen GDH** no cross reactivity against organism, pathogens, substances that could cause infections:

Cross reactivity										
Virus	Adenovirus	-	Coronavirus	-	Influenza B	-	Norovirus GII	-	Rotavirus	-
	Astrovirus	-	Influenza A	-	Norovirus GI	-	Respiratory Syncytial Virus	-		
Bacteria	Campylobacter coli	-	Clostridium Haemolyticum	-	Helicobacter pylori	-	Salmonella typhimurium	-	Streptococcus pyogenes	-
	Campylobacter jejuni	-	Clostridium Novyi	-	Legionella pneumophila	-	Shigella boydii	-	Yersinia enterocolitica O:3	-
	Clostridium difficile Toxin A	-	Clostridium Tetani	-	Listeria monocytogenes	-	Shigella dysenteriae	-	Yersinia enterocolitica O:9	-
	Clostridium difficile Toxin B	-	Clostridium Septicum	-	Peptostreptococcus anaerobius	-	Shigella flexneri	-		
	Clostridium perfringens	-	Escherichia coli O:157	-	Salmonella enteritidis	-	Shigella sonnei	-		
	Clostridium bifementans	-	Escherichia coli O:026	-	Salmonella paratyphi	-	Staphylococcus aureus	-		
	Clostridium Butyricum	-	Escherichia coli O:111	-	Salmonella typhi	-	Streptococcus pneumococcal	-		
Parasites	Cryptosporidium Parvum	-	Entamoeba histolytica	-	Giardia lamblia	-				
Others	Pig haemoglobin	-	Bovine haemoglobin		Bovine Lactoferrin	-	Bovine Transferrin	-		

Table 5 Cross reactivity results.

A specificity assay was performed for **CerTest Clostridium difficile antigen GDH**. This concluded that the test, could detects the following antigens: *C. sporogenes* (CECT 485) and *C. botulinium* (CECT 551).

In this specificity assay, *Clostridium sordellii* (ATCC 9714) was evaluated and no cross reactivity was found. In few intercomparison studies **CerTest Clostridium difficile antigen GDH** have a positive signal with this pathogen, but the nature of this positive result could not be established.

E. Interferences

An evaluation was performed to determine the possible interferences of **CerTest Clostridium difficile antigen GDH**; no interferences against the substances tested were detected:

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metronidazole	3.0 mg/mL	No interference	Loperamide hydrochloride (Fortasec)	0.15 mg/mL	No interference
	Ampicillin	3.0 mg/mL	No interference	Heparin (Hibar)	350.0 IU/mL	No interference
	Oseltamivir	0.3 mg/mL	No interference	Almagato (Almax)	3.0 mg/mL	No interference
	Amantadine	0.3 mg/mL	No interference	Fosfamycin (Monurol)	$3.0 \cdot 10^{-3}$ mg/mL	No interference
	Ribavirin	3.0 mg/mL	No interference	Acetylcysteine (Fluimucil)	3.0 mg/mL	No interference
	Codeine (Toseina)	0.2 mg/mL	No interference	Dexketoprofen trometamol (Enantyum)	0.3 mg/mL	No interference
	Benzocaine (Angileptol)	0.3 mg/mL	No interference	Levofloxacin	3.0 mg/mL	No interference
	Cloperastine (Flutox)	0.3 mg/mL	No interference	Ciprofloxacin	0.3 mg/mL	No interference
	Carbocisteine mucolítico (Iniston)	3.0 mg/mL	No interference	Rifampicin (Rifaldin)	0.3 mg/mL	No interference
	Loratadine	0.3 mg/mL	No interference	Phenoxymethylpenicillin potassium	3.0 mg/mL	No interference
	Dexchlorpheniramine (Polaramine)	0.3 mg/mL	No interference	Ambroxol hydrochloride (Mucosan)	0.3 mg/mL	No interference
	Ebastine (Ebastel)	3.0 mg/mL	No interference	Macrogol 3350 (Movicol)	3.0 mg/mL	No interference
	Acetyl Salicylic (Adiro)	0.3 mg/mL	No interference	Lysine Carbocysteinatate (Pectox)	$3.0 \cdot 10^{-2}$ mg/mL	No interference
	Ibuprofen (Espidifen)	3.0 mg/mL	No interference	Hydroxyzine dihydrochloride	0.3 mg/mL	No interference
	Paracetamol (Dolocartil)	5.0 mg/mL	No interference	Lorazepam	$3.0 \cdot 10^{-3}$ mg/mL	No interference
	Metamizole (Nolofil)	5.0 mg/mL	No interference	Amoxicillin	3.0 mg/mL	No interference
	Prednisone	0.3 mg/mL	No interference	Mercaptopurine	0.3 mg/mL	No interference
	Omeprazole	$2.0 \cdot 10^{-3}$ mg/mL	No interference	Biotine	100.0 µg/mL	No interference
	Naso GEL	0.9 mg/mL	No interference	Sore Throat Phenol spray	0.5 mg/mL	No interference
	CVS Nasal Spray (Cromolyn)	4 mg/mL	No interference	Tobramycin	0.3 mg/mL	No interference
Afrin (Oxymetazoline)	$5.0 \cdot 10^{-2}$ mg/mL	No interference	Mupirocin	$2.5 \cdot 10^{-2}$ mg/mL	No interference	
CVS Nasal Drops (Phenylephrine)	10.0 mg/mL	No interference	Fluticasone Propionate	$5.0 \cdot 10^{-2}$ mg/mL	No interference	
ZICAM	1.0 mg/mL	No interference	Homeopathic	DIL 1/10	No interference	

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
	Phenylpropanolamine	0.3 mg/mL	No interference			
Endogenous	Human haemoglobin	5.0 µg/mL	No interference	Human lactoferrin	5.0 µg/mL	No interference
	Human transferrin	0.5 µg/mL	No interference	Mucine	0.5 %	No interference
	Human calprotectin	5.0 µg/mL	No interference	Human blood	5.0 %	No interference

Table 6 Interferences results.

F. Repeatability and Reproducibility

Following the study of repeatability and reproducibility of of **CerTest Clostridium difficile antigen GDH**, performed using different internal samples, negative and positive, no differences have been observed within the evaluations.

For information purposes only

ANNEX 1

KIT 20 CARD *CLOSTRIDIUM DIFFICILE* ANTIGEN GDH. Catalogue reference: GD820001V

Annex for the following reference:

PRODUCT	REFERENCE
CerTest <i>Clostridium difficile</i> antigen GDH Card (20 tests kit with stool vials)	GD820001V

Table A 1. Reference for kit 20 card of CerTest *Clostridium difficile* antigen GDH.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
-20x CerTest <i>Clostridium difficile</i> antigen GDH (card test format) -20x Stool collection tubes with diluent - 1x IFU Instructions for use (including Quick guide use).	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Specimen collection container


Table A 1.1 Materials necessities for test performance.

II. Test procedures (samples)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the **stool collection tube** in order to assure good sample dispersion.
2. Remove **CerTest *Clostridium difficile* antigen GDH** (card test) from its sealed bag just before using it.
3. Take the **stool collection tube**, cut the end of the cap (1) and dispense 3 drops in the circular window marked with the letter S (2). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

III. Interpretation of results

(please, refer to the illustration below)

	CerTest <i>Clostridium difficile</i> antigen GDH	Interpretation of the results
1.	NEGATIVE	<p>There is no <i>Clostridium difficile</i> presence in the sample. No infection caused by <i>Clostridium difficile</i>. A negative result should be followed up with additional laboratory techniques to confirm the results.</p>
	GREEN	
2.	POSITIVE	<p>There is <i>Clostridium difficile</i> antigen presence in the sample. Possible infection caused by <i>Clostridium difficile</i>. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.</p>
	GREEN-RED	
3.	INVALID	<p>Invalid result, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.</p>
	Any other result	

Table A 1.2 CerTest *Clostridium difficile* antigen GDH performance interpretation of the results.

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control

INTERNAL CONTROL: control line

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

For information purposes only

ANNEX 2

KIT 20 CARD *CLOSTRIDIUM DIFFICILE* ANTIGEN GDH. Catalogue reference: GD820001VC

Annex for the following reference:

PRODUCT	REFERENCE
CerTest <i>Clostridium difficile</i> antigen GDH Card (20 tests kit with stool vials + PC)	GD820001VC

Table A 2. Reference for kit 20 card of CerTest *Clostridium difficile* antigen GDH with positive control.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
-20x CerTest <i>Clostridium difficile</i> antigen GDH (card test format) -20x Stool collection tubes with diluent. -1x CerTest <i>C. difficile</i> antigen GDH Positive Control -1x Testing tube -1x Disposable pipette -1x IFU Instructions for use (including Quick guide use).	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Specimen collection container


Table A 2.1 Materials necessities for test performance.

II. Test procedures (samples)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the **stool collection tube** in order to assure good sample dispersion.
2. Remove **CerTest *Clostridium difficile* antigen GDH** (card test) from its sealed bag just before using it.
3. Take the **stool collection tube**, cut the end of the cap (1) and dispense 3 drops in the circular window marked with the letter S (2). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

III. Interpretation of results

(please, refer to the illustration below)

CerTest <i>Clostridium difficile</i> antigen GDH		Interpretation of the results	
1.	NEGATIVE		<p>There is no <i>Clostridium difficile</i> presence in the sample. No infection caused by <i>Clostridium difficile</i>.</p> <p>A negative result should be followed up with additional laboratory techniques to confirm the results.</p>
	GREEN		
2.	POSITIVE		<p>There is <i>Clostridium difficile</i> antigen presence in the sample. Possible infection caused by <i>Clostridium difficile</i>.</p> <p>The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.</p> <p>A positive result should be followed up with additional laboratory techniques to confirm the results.</p>
	GREEN-RED		
3.	INVALID		<p>Invalid result, we recommend repeating the assay using the same sample with another test.</p> <p>Any result read after more than 10 minutes will be considered invalid.</p>
	Any other result		

Table A.2.2 CerTest *Clostridium difficile* antigen GDH performance interpretation of the results.

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control

CerTest Clostridium difficile antigen GDH (card test) presents two controls:

INTERNAL CONTROL: control line

Internal control: internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROL: POSITIVE CONTROL SWAB

External control: CerTest C. difficile antigen GDH Positive Control (positive control swab).

A. External control. CerTest C. difficile antigen GDH Positive Control

CerTest C. difficile antigen GDH Positive Control is an external quality control for **CerTest Clostridium difficile antigen GDH**. Use **CerTest C. difficile antigen GDH Positive Control** to check that the extraction reagents and the test are working properly.


Store as packaged in the sealed pouch at 2-30°C. **CerTest C. difficile antigen GDH Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest C. difficile antigen GDH Positive Control** must remain in the sealed pouch until use. Do not freeze.

A.1 Positive control procedure

Allow tests and CerTest C. difficile antigen GDH Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Add 15 drops (1) of diluent from the **stool collection tube** into the **testing tube** and immediately introduce **CerTest C. difficile antigen GDH Positive Control**.
2. Mix the solution by rotating **CerTest C. difficile antigen GDH Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from **CerTest C. difficile antigen GDH Positive Control**, squeezing the sides of the **testing tube** as **CerTest C. difficile antigen GDH Positive Control** is withdrawn (2). Discard **CerTest C. difficile antigen GDH Positive Control**.

3. Remove **CerTest Clostridium difficile antigen GDH** (card test) from its sealed bag just before using it.
4. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter S (3).
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Positive control could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1.1 Interpretation positive control procedure

CerTest <i>C. difficile</i> antigen GDH Positive Control		Interpretation of the results
1.	POSITIVE	 <p>A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window during the test performance due to presence of GDH antigen in CerTest <i>C. difficile</i> antigen GDH Positive Control.</p>
	GREEN-RED	
2.	INVALID	<p>Total absence of the control coloured line (GREEN), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new CerTest <i>C. difficile</i> antigen GDH Positive Control.</p> <p>If the situation still persists, discontinue using the test kit and contact your local distributor.</p> <p>Any result read after more than 10 minutes will be considered invalid.</p>
	Any other result	

Table A 2.3 CerTest *C. difficile* antigen GDH Positive Control performance interpretation of the results.

ANNEX 3

KIT 20 CARD *CLOSTRIDIUM DIFFICILE* ANTIGEN GDH. Catalogue reference: GD820001VCC

Annex for the following reference:

PRODUCT	REFERENCE
CerTest <i>Clostridium difficile</i> antigen GDH Card (20 tests kit with stool vials + PC + NC)	GD820001VCC

Table A 3. Reference for Kit 20 card of CerTest *Clostridium difficile* GDH antigen with positive and negative control.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
-20x CerTest <i>Clostridium difficile</i> antigen GDH (card test format) -20x Stool collection tubes with diluent. -1x CerTest <i>C. difficile</i> antigen GDH Positive Control -1x CerTest Negative Control -2x Testing tube -2x Disposable pipette -1x IFU Instructions for use (including Quick guide use).	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Specimen collection container

Table A 3.1 Materials necessities for test performance.


II. Test procedures (samples)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the **stool collection tube** in order to assure good sample dispersion.
2. Remove **CerTest *Clostridium difficile* antigen GDH** (card test) from its sealed bag just before using it.
3. Take the **stool collection tube**, cut the end of the cap (1) and dispense 3 drops in the circular window marked with the letter S (2). Avoid adding solid particles with the liquid.

4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

Wait 10 minutes

III. Interpretation of results

(please, refer to the illustration below)

	CerTest <i>Clostridium difficile</i> antigen GDH	Interpretation of the results
1.	NEGATIVE	<p>There is no <i>Clostridium difficile</i> presence in the sample. No infection caused by <i>Clostridium difficile</i>.</p> <p>A negative result should be followed up with additional laboratory techniques to confirm the results.</p>
	GREEN	
2.	POSITIVE	<p>There is <i>Clostridium difficile</i> antigen presence in the sample. Possible infection caused by <i>Clostridium difficile</i>.</p> <p>The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.</p> <p>A positive result should be followed up with additional laboratory techniques to confirm the results.</p>
	GREEN-RED	
3.	INVALID	<p>Invalid result, we recommend repeating the assay using the same sample with another test.</p> <p>Any result read after more than 10 minutes will be considered invalid.</p>
	Any other result	

Table A 3.2 CerTest *Clostridium difficile* antigen GDH performance interpretation of the results.

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control

CerTest *Clostridium difficile* antigen GDH (card test) presents three controls:

INTERNAL CONTROL: control line

Internal control: internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROL: POSITIVE CONTROL AND NEGATIVE CONTROL SWAB:

External positive control: CerTest *C. difficile* antigen GDH Positive Control (positive control swab).

External negative control: CerTest Negative Control (negative control swab).

A. External positive control. CerTest *C. difficile* antigen GDH Positive Control

CerTest *C. difficile* antigen GDH Positive Control is an external quality control for **CerTest *Clostridium difficile* antigen GDH**. Use **CerTest *C. difficile* antigen GDH Positive Control** to check that the extraction reagents and the test are working properly.


Store as packaged in the sealed pouch at 2-30°C, **CerTest *C. difficile* antigen GDH Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest *C. difficile* antigen GDH Positive Control** must remain in the sealed pouch until use. Do not freeze.

A.1 Positive control procedure

Allow tests and CerTest *C. difficile* antigen GDH Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Add 15 drops (1) of diluent from the **stool collection tube** into the **testing tube** and immediately introduce **CerTest *C. difficile* antigen GDH Positive Control**.

- Mix the solution by rotating **CerTest C. difficile antigen GDH Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from **CerTest C. difficile antigen GDH Positive Control**, squeezing the sides of the **testing tube** as **CerTest C. difficile antigen GDH Positive Control** is withdrawn (2). Discard **CerTest C. difficile antigen GDH Positive Control**.
- Remove **CerTest Clostridium difficile antigen GDH** (card test) from its sealed bag just before using it.
- Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the circular window marked with the letter S (3).
- Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Positive control could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1.1 Interpretation positive control procedure

	CerTest C. difficile antigen GDH Positive Control	Interpretation of the results
1.	<p>POSITIVE</p> <p>GREEN-RED</p>	A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window during the test performance due to presence of GDH antigen in CerTest C. difficile antigen GDH Positive Control .
2.	<p>INVALID</p> <p>Any other result</p>	Total absence of the control coloured line (GREEN), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the

CerTest <i>C. difficile</i> antigen GDH Positive Control		Interpretation of the results
		<p>control with a new test and a new CerTest <i>C. difficile</i> antigen GDH Positive Control.</p> <p>If the situation still persists, discontinue using the test kit and contact your local distributor.</p> <p>Any result read after more than 10 minutes will be considered invalid.</p>

Table A 3.3 CerTest *C. difficile* antigen GDH Positive Control performance interpretation of the results.


B. External negative control. CerTest Negative Control

CerTest Negative Control is an external quality control for **CerTest *Clostridium difficile* antigen GDH** (card test). Use **CerTest Negative Control** to check that the extraction reagents and the test are working properly. Store as packaged in the sealed pouch at 2-30°C, **CerTest Negative Control** is stable through the expiration date printed on the sealed pouch. **CerTest Negative Control** must remain in the sealed pouch until use. Do not freeze.

B.1 Negative control procedure

Allow tests and CerTest Negative Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Add 15 drops (1) of diluent from the **stool collection tube** into the **testing tube** and immediately put **CerTest Negative Control**.
2. Mix the solution by rotating **CerTest Negative Control** forcefully against the side of the **testing tube** at least 15 seconds. Extract as much liquid as possible from **CerTest Negative Control**, squeezing the sides of the **testing tube** as **CerTest Negative Control** is withdrawn (2). Discard **CerTest Negative Control**.
3. Remove **CerTest *Clostridium difficile* antigen GDH** (card test) from its sealed bag just before using it.
4. Use a separate **disposable pipette** and test for each control. Dispense 3 drops from the **testing tube** in the corresponding circular window (3).
5. **Read the result at 10 minutes**. Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Negative controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

B.1.1. Interpretation negative control procedure

CerTest Negative Control		Interpretation of the results
1.	VALID RESULT	 <p>Only A GREEN line (control line (C)) appear across the results window during the test performance.</p>
	GREEN (only control band appear)	
2.	INVALID	 <p>Total absence of the control coloured line (GREEN), regardless the appearance or not of the test line (RED), or present of both lines (RED test and GREEN control lines). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new CerTest Negative Control.</p> <p>If the situation still persists, discontinue using the test kit and contact your local distributor.</p> <p>Any result read after more than 10 minutes will be considered invalid</p>
	Any other result	

Table A 3.4 CerTest Negative Control performance interpretation of the results.

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11. Reference/Bibliografía

1. D. Rodríguez-Pardo et al. **Infecciones producidas por Clostridium difficile**. *Enfermedades Infecciosas y Microbiología Clínica*, 2013, 31(4): 254-263.
2. Fitzsimmons. **Independent evaluation of a GDH enzyme immunoassay for C. difficile screening**. *Biomed Sci.* 2011;(October):698–700.
3. Vaishnavi C. **Clinical spectrum & pathogenesis of Clostridium difficile associated diseases**. *Indian J Med Res.* 2010;131(4):487–99.
4. Wren MWD, Sivapalan M, Kinson R, Shetty NP. **An evaluation of tests for faecal toxin , glutamate dehydrogenase , lactoferrin and toxigenic culture in the diagnostic laboratory**. *Br J Biomed Sci.* 2009;66(1):1–5.
5. Ofori E, Ramai D, Dhawan M, Mustafa F, Gasperino J, Reddy M. **Community-acquired Clostridium difficile: epidemiology, ribotype, risk factors, hospital and intensive care unit outcomes, and current and emerging therapies**. *J Hosp Infect.* 2018;99(4):436–42.
6. Guh AY, Mu Y, Winston LG, Johnston H, Olson D, Farley MM, et al. **Trends in U.S. Burden of Clostridioides difficile Infection and Outcomes**. *N Engl J Med.* 2021;382(14):1320–30.
7. Services U. D of H and H. **Antibiotic resistance threats in the United States**. *Centers Dis Control Prev.* 2019;1–113.

12. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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