Calprotectin TURBILATEX[®] Combo

Calprotectin latex turbidimetric assay

Cat. TL-022CP100ED, TL-022CP200ED, TL-022CP400ED, TL-022CP800ED

IU-TL – 022CP100ED-200ED-400ED-800ED en-es v.08

CONTENT/CONTENIDO

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ENGLISH

INTENDED USE

Calprotectin Turbilatex® is a latex turbidimetric assay **only for the quantitative detection of calprotectin in human stool samples** (not to be used for body fluid as blood, serum, plasma, urine, cerebrospinal fluid, oral fluid, synovial fluid or empyema fluid).

This assay is simple and widely applicable. This product is optimized for several automated analyser.

For professional *in vitro* diagnostic use only.

SYNTHESIS

Calprotectin is a surrogate marker of neutrophil turnover and is elevated in various inflammatory conditions. It has been studied in numerous infectious and inflammatory diseases, especially in gastroenterology, where it is elevated in several gastrointestinal conditions. Calprotectin is detectable in stool and is highly resistant to proteolytic degradation in the presence of calcium, making its measurement in stool samples useful for providing diagnostic and prognostic information for various common organic diseases (1). Studies have shown that faecal calprotectin levels correlate with colorectal and intestinal inflammation (2), and it is considered a useful marker for inflammatory bowel disease (IBD) (3).

IBD refers to a group of gastrointestinal inflammatory diseases, with Crohn's disease (CD) and ulcerative colitis (UC) being the main types (4). As these are inflammatory disorders, a neutrophil-derived protein like calprotectin can be valuable for diagnosis and monitoring (1). Both CD and UC are characterized by periods of remission and episodes of relapse, often caused by acute intestinal inflammation (5), leading to symptoms like diarrhoea, abdominal pain, and rectal bleeding (4). However, most patients with such symptoms have irritable bowel syndrome (IBS) rather than IBD, with IBS showing normal colonoscopy results (4). Calprotectin detection can reduce the need for invasive colonoscopies (4).

PRINCIPLE

Calprotectin Turbilatex® Combo is a latex turbidimetric assay **only for the quantitative detection of Calprotectin in human stool samples**. This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use of the test is **exclusively to differentiate IBD patients with inflammation from IBD patients without inflammation and from irritable bowel syndrome (IBS)**.

Calprotectin latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific

reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample, it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at determinate wavelength and is proportional to the quantity of antigen contained in the sample.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in Turbidimetric technique, autoanalyzer use and *in vitro* diagnostic procedures.
- The analyser should be ready before performing any assay.
- Read the instructions for use carefully before use the kit.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Do not freeze to maintain the integrity and performance of the test.
- Do not mix screw caps from different reagent vials.
- Components provided in the kit are approved for use in **Calprotectin Turbilatex® Combo**. Do not mix or use with components from other product kit or lots. Do not use with other commercial kit components (from another supplier). The test's satisfactory performance is ensured only when components from the same lot of **Calprotectin Turbilatex® Combo** are used.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples.
- **Calprotectin Turbilatex® Combo** do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which

are in concentrations higher than the value established in the mentioned regulation for their declaration. Material Safety Data Sheet are not included with this device, nevertheless it could be consulted upon request to Certest Biotec.

- **Calprotectin Turbilatex® Combo** has been validated with the analyser Biolis from Tokyo Boeki. For different analysers, Calprotectin should be measured according to the specific application parameters for each specific chemistry analyzer. More information about the application notes available can be found on our website.
- Homogenise the reagents before use, make sure that the reagents are free of bubbles.
- Evaporation of calibrators and controls can lead to incorrect results. Follow the storage recommendations detailed in the following section.
- If the values obtained for the controls are not valid, repeat the quality control measurement with new controls. If the control values are still invalid, recalibrate the equipment.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- The solutions should be discarded in a proper container after testing following local regulations.

STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). **Do not freeze.**

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight and conserved at 2-8°C.

REAGENTS

Materials provided

<u>References</u>	<u>Cat.:</u> <u>TL-022CP100ED</u>	<u>Cat.:</u> <u>TL-022CP200ED</u>	<u>Cat.:</u> <u>TL-022CP400ED</u>	<u>Cat.:</u> <u>TL-022CP800ED</u>
- TL-022CP01 Calprotectin Turbilatex® Reagent 1	1x27mL	2 x 27mL	4 x 27mL	8 x 27mL
- TL-022CP02 Calprotectin Turbilatex® Reagent 2	1 x 5mL	1 x 8mL	2 x 8mL	4 x 8mL
- TL-022CP70/ TL-022CP71/ TL-022CP72- TL-022CP73/ TL-022CP74/ TL-022CP75 Calprotectin Turbilatex® Calibrators		6 x 1mL		
- TL-022CP08/ TL-022CP09 Calprotectin Turbilatex® Control 1&2		4 x 1mL		

Materials and instruments that could be required but no provided

- Specimen collection container.
- Disposable gloves and laboratory equipment.

- Automated analyser.
- MST-0019U Universal Turbilatex® Sample Collection Vial.

SPECIMEN COLLECTION AND PREPARATION

Collect enough quantity of human solid stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. If not immediately tested, freeze the stored samples at -20 °C maximum 6 months. In this case, the sample will be totally thawed and brought to room temperature (15-30°C) before testing. Homogenize stool samples as thoroughly as possible prior to preparation.

To process the collected stool samples using catalogue reference MST-0019U follow the indication of the instruction for use of the product.

ANALYTICAL PROCEDURE

R1 and R2 are ready to use.

Preparation of the calibration curve


For calibration use only the following materials:

TL-022CP70, TL-022CP71, TL-022CP72, TL-022CP73, TL-022CP74 and TL-022CP75, Calprotectin Turbilatex® Calibrator vials: liquid calibrator containing recombinant human calprotectin at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a month is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

Calprotectin Turbilatex® Calibrator vials (calibrator 0 to 5) are ready to use.

Picture 1. Prepared calibrators for performing the calibration curve (5 points and blank)

Quality control

For quality control only use the following materials:


TL-022CP08 Calprotectin Turbilatex® Control 1 and TL-022CP09 Calprotectin Turbilatex® Control 2: liquid controls at two different concentrations of recombinant human calprotectin. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.

Calprotectin Turbilatex® Control vials (Control 1 and Control 2) are ready to use.

Control vials

Picture 2. Prepared controls C1 and C2 (check control concentration on the label)

INTERPRETATION OF RESULTS

Positive results: higher or equal than the cut-off fixed by the clinical lab.

Recommended: 50 µg of hCp/g of stool for diagnostic procedures and 200 µg of hCp/g of stool for screening procedures.

Positive results determine the abnormal presence of human Calprotectin (hCp) in stool samples.

Calprotectin Turbilatex® Combo cut-off value:

Calprotectin concentration values lower than 50 µg of hCp/g of stool are considered normal values and that is not indicative of an inflammation of gastrointestinal tract.

Calprotectin concentration values between 50 and 200 µg hCp/g of stool are considered as an abnormal presence being indicative of mild inflammation of gastrointestinal tract, therefore monitoring and follow-up of the patient is recommended.

Calprotectin concentration values higher than 200 µg of hCp/g of stool are indicative of a severe inflammation of gastrointestinal tract.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Limit of detection (LOD): 7 µg hCp/g of stool. The lower limit of detection of **Calprotectin Turbilatex® Combo** was determined on 20 samples and 2 sample replicates as the mean value +2· SD.

Limit of quantification (LOQ): 20 µg hCp/g of stool. The lower limit of quantification is defined as the lowest actual amount of analysis that can be reliably detected; imprecision is < 20% as CV%.

Measuring range

Calprotectin Turbilatex® Combo assay measuring range is 20-3125 µg hCp/g of stool. Samples higher concentrated than 1500 µg hCp/g of stool must be diluted for proper quantification by the user, using additional sample buffer.

Prozone

Studies have been made up to a concentration of 25000 µg of hCp/g of stool and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Samples with concentrations up to 3125 µg of hCp/g can be measured without inhibitory prozone effect.

Clinical sensitivity and specificity

The following tables provide information about the evaluations **Calprotectin Turbilatex® Combo**, has undergone with a cut-off of 50 µg of hCp/g of stool:

N°	Evaluated test	Sample type
1	Calprest ELISA, Eurospital	Human stool samples
2	Calprest ELISA, Eurospital	Human stool samples
3	BÜHLMANN fCAL® turbo, Bühlmann Diagnostics	Human stool samples
4	BÜHLMANN fCAL® turbo, Bühlmann Diagnostics	Human stool samples

Table 1 Resume of Clinical evaluations. Test and sample type.

N°	TP	FP	FN	TN	N	Sens. (%) CI (%)	Spec. (%) CI (%)	PPV (%) CI (%)	NPV (%) CI (%)	LR+ CI	LR- CI
1	16	0	1	38	55	94.1 (71.3-99.9)	100.0 (79.4-100.0)	100.0 (90.7-100.0)	97.4 (96.5-99.9)	71.50 (4.58-1127)	0.084 (0.018-0.391)
2	42	6	0	52	100	100.0 (91.6-100.0)	89.7 (78.8-96.1)	87.5 (74.8-95.3)	100.0 (93.2-100.0)	8.97 (4.341-18.54)	0.013 (0.001-0.206)
3	51	7	3	39	100	94.4 (84.6-98.8)	84.8 (71.1-93.7)	87.9 (76.7-95.0)	92.9 (80.5-98.5)	6.206 (3.128-12.31)	0.066 (0.022-0.198)
4	42	5	1	62	110	97.7 (87.7-99.9)	92.5 (83.4-97.5)	89.4 (76.9-96.5)	98.4 (91.5-100.0)	13.09 (5.63-30.45)	0.025 (0.004-0.175)
Total	151	18	5	191	365	96.8 (92.7-99.0)	91.4 (86.7-94.8)	89.3 (83.7-93.6)	97.4 (94.1-99.2)	11.24 (7.22-17.49)	0.035 (0.015-0.083)

Table 2 Clinical evaluations and global results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI), Likelihood positive ratio (LR+), Likelihood negative ratio (LR-) for *Calprotectin Turbilatex® Combo*.

The results showed a high sensitivity and specificity to detect human calprotectin (hCp) using *Calprotectin Turbilatex® Combo*.

Interferences

An evaluation was performed to determine the interferences of *Calprotectin Turbilatex® Combo*. No interferences were founded against other substances tested: BSA, Hemin, Ascorbic acid.

Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against other faecal markers occasionally present in faeces, such as: bovine and pig haemoglobin, bovine transferrin, bovine lactoferrin and human haemoglobin, transferrin and lactoferrin.

Precision

An evaluation was performed to determine the accuracy of *Calprotectin Turbilatex® Combo*. *Calprotectin Turbilatex® Combo* shows high in its reproducibility results on known concentration of calprotectin.

Linearity

Calprotectin Turbilatex® Combo using calibrator kit is linear in the calibration range of 36-1500 µg hCp/g of stool.

The analytical performance of *Calprotectin Turbilatex® Combo* has been validated using the Tokyo Boeki Biolis 24i analyser. For more information about other compatible equipment,

please refer to the corresponding application. The different turbidimetric equipment with which our turbidimetric products are compatible can be found in **Annex I**.

LIMITATIONS

1. **Calprotectin Turbilatex® Combo** should **only be used in human stool samples (not to be used for body fluid such as blood, serum, plasma, urine, cerebrospinal fluid, oral fluid, synovial fluid or empyema fluid)**. The use of other samples has not been established. Proper faecal specimens must be obtained.
2. Positive results determine the presence of hCp in faecal samples. A positive result should be followed up with additional diagnostic invasive procedures, a colonoscopy and a biopsy in order to confirm the diagnosis and to establish the inflammation extent.
3. If symptoms or situation still persist, calprotectin determination should be carried out invasive techniques. Negative results do not exclude IBD with inflammation, some diseases such as celiac sprue and microscopic colitis polyps that mainly involve mononuclear inflammation (6,7).
4. Stool samples from patients with non-steroidal anti-inflammatory drug (NSAID) treatment could show positive result (7).
5. Neonatal faecal calprotectin levels have been reported higher than those in normal children with a mean of 167 µg/g (range 22-860 µg/g) (8).
6. Patients with agranulocytosis may produce false-negative results due to low neutrophil production in the bone marrow, which affects calprotectin detection (9).
7. Patients with Inflammatory Bowel Disease (IBD) transition between active (inflammatory) and inactive phases of the disease. It is important to consider these disease activity changes when interpreting **Calprotectin Turbilatex® Combo** assay results (9).
8. Other intestinal disorders, such as gastrointestinal infections and colorectal cancer, can also elevate calprotectin levels. Therefore, a positive result on the **Calprotectin Turbilatex® Combo** assay cannot be used as the sole criterion for diagnosing active IBD (6).

EXPECTED VALUES

Inflammatory bowel diseases (IBD), including Crohn's disease and ulcerative colitis, are chronic conditions affecting millions worldwide. The prevalence varies significantly by region, with rates in North America, Oceania and Europe reaching up to 0.3% of the population, which makes IBD one of the most common gastrointestinal disorders in these areas. The rise in cases

is attributed partly to greater awareness, environmental factors, and potentially genetic predispositions in certain populations (10).

Some studies established equal or higher 50 µg hCp/g faeces as cut-off value to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

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1. Poullis A, Foster R, Mendall MA, Fagerhol MK. **Emerging role of calprotectin in gastroenterology. Vol. 18, Journal of Gastroenterology and Hepatology (Australia).** Blackwell Publishing; 2003. p. 756–62.
2. Angriman I, Scarpa M, D'Inca R, Basso D, Ruffolo C, Polese L, et al. **Enzymes in feces: Useful markers of chronic inflammatory bowel disease.** Clinica Chimica Acta. 2007 May 2;381(1):63–8.
3. Mayer EA. **Functional Gastrointestinal Disorders: Irritable Bowel Syndrome, Dyspepsia, and Functional Chest Pain of Presumed Esophageal Origin.** Goldman's Cecil Medicine: Twenty Fourth Edition. 2012;1:868–74.
4. Khaki-Khatibi F, Qujeq D, Kashifard M, Moein S, Maniati M, Vaghari-Tabari M. **Calprotectin in inflammatory bowel disease.** Vol. 510, Clinica Chimica Acta. Elsevier B.V.; 2020. p. 556–65.
5. Tibble JA, Sigthorsson G, Bridger S, Fagerhol MK, Bjarnason I. **Surrogate markers of intestinal inflammation are predictive of relapse in patients with inflammatory bowel disease.** Gastroenterology. 2000;119(1):15–22.
6. Pathirana WGW, Chubb SP, Gillett MJ, Vasikaran SD. **Faecal Calprotectin.** Clin Biochem Rev. 2018;39(3):77.
7. Klingberg E, Strid H, Ståhl A, Deminger A, Carlsten H, Öhman L, et al. **A longitudinal study of fecal calprotectin and the development of inflammatory bowel disease in ankylosing spondylitis.** Arthritis Res Ther. 2017 Feb 2;19(1):1–9.
8. Campeotto F, Butel MJ, Kalach N, Derrieux S, Aubert-Jacquín C, Barbot L, et al. **High faecal calprotectin concentrations in newborn infants.** Arch Dis Child Fetal Neonatal Ed. 2004 Jul;89(4).
9. Fukunaga S, Kuwaki K, Mitsuyama K, Takedatsu H, Yoshioka S, Yamasaki H, et al. **Detection of calprotectin in inflammatory bowel disease: Fecal and serum levels and immunohistochemical localization.** Int J Mol Med. 2018;41(1):107–18.
10. Ng SC, Shi HY, Hamidi N, Underwood FE, Tang W, Benchimol EI, et al. **Worldwide incidence and prevalence of inflammatory bowel disease in the 21st century: a systematic review of population-based studies.** The Lancet. 2017 Dec 23;390(10114):2769–78.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO *IN VITRO*

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website: www.certest.es

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