Pancreatic Elastase TURBILATEX[®] Combo

Pancreatic Elastase latex turbidimetric assay

Cat. TL-022EL100ED, TL-022EL200ED, TL-022EL400ED, TL-022EL800ED

IU-TL - 022EL100ED-200ED-400ED-800ED en-es v.06

CONTENT/CONTENIDO

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ENGLISH

INTENDED PURPOSE

Pancreatic Elastase Turbilatex® is a latex turbidimetric assay **only for the quantitative detection of Pancreatic Elastase E1 in human stool samples.**

This assay is simple and widely applicable. This product is optimized for several automated analyser.

Test results should **exclusively be used to evaluate exocrine pancreatic function in stool samples.**

For professional *in vitro* diagnostic use only.

SYNTHESIS

Pancreatic elastase Turbilatex® Combo is a non-invasive *in vitro* diagnostic method, based on monoclonal antibodies and used to monitor exocrine pancreatic function.

Exocrine pancreatic insufficiency (EPI) is a condition characterized by the inability of the pancreas to produce digestive enzymes in the required amount, leading to poor digestion and malabsorption of nutrients. This can be caused by various conditions, including chronic pancreatitis, pancreatic cancer, and other diseases that affect the structure and function of the pancreas (1,2). A significant percentage of patients with EPI may experience symptoms during childhood, and many may develop this condition over their lifetime (3).

EPI is defined as a reduction in the activity of pancreatic enzymes in the intestinal lumen below the threshold necessary for normal digestion (4). Patients with EPI often require pancreatic enzyme replacement therapy (PERT), which involves daily treatment and can be costly (5). Continuous monitoring is essential, as EPI can develop without evident symptoms, complicating diagnosis (2).

Pancreatic elastase is an enzyme produced by the pancreas and is considered an effective non-invasive biomarker for evaluating exocrine pancreatic function. Pancreatic elastase is stable in stool samples and does not degrade during intestinal transit, making its measurement straightforward. Classic symptoms of EPI include weight loss, gas, bloating, and steatorrhea (6). However, many patients may be asymptomatic, highlighting the need for accurate diagnostic methods, as available tests may only yield positive results in advanced stages of insufficiency.

PRINCIPLE

Pancreatic Elastase Turbilatex® is a latex turbidimetric assay **only for the quantitative detection of Pancreatic elastase E1 in human stool samples**. This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use of the test is **exclusively be used to evaluate Pancreatic elastase levels in stool samples associated to Pancreatic insufficiency**.

Pancreatic Elastase latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and monoclonal antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample, it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at the determined wavelength and is proportional to the quantity of antigen contained in the sample.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in Turbidimetric technique, autoanalyzer use and *in vitro* diagnostic procedures.
- The analyser should be ready before performing any assay.
- Read the instructions for use carefully before use the kit. Do not perform the test without fully understanding the information on procedures, safety precautions and limitations described in this instruction for use.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Do not freeze to maintain the integrity and performance of the test.
- Do not mix screw caps from different reagent vials.
- The kit does not contain materials of human origin.
- Components provided in the kit are approved for use in **Pancreatic Elastase Turbilatex® Combo**. Do not mix or use with components from other product kit or lots. Do not use with other commercial kit components (from another supplier). The test's satisfactory

performance is ensured only when components from the same lot of **Pancreatic Elastase Turbilatex® Combo** are used.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples.
- **Pancreatic Elastase Turbilatex® Combo** do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration. Material Safety Data Sheet are not included with this device, nevertheless it could be consulted upon request to Certest Biotec.
- **Pancreatic Elastase Turbilatex® Combo** has been validated with the analyser Biolis from Tokyo Boeki. For different analysers, Pancreatic elastase should be measured according to the specific application parameters for each specific chemistry analyzer. More information about the application notes available can be found on our website.
- Homogenise the reagents before use, make sure that the reagents are free of bubbles.
- Evaporation of calibrators and controls can lead to incorrect results. Follow the storage recommendations detailed in the following section.
- If the values obtained for the controls are not valid, repeat the quality control measurement with new controls. If the control values are still invalid, recalibrate the equipment.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- The solutions should be discarded in a proper container after testing following local regulations.

STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). **Do not freeze.**

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight and conserved at 2-8°C.

REAGENTS

Materials provided

References	Cat.: TL-022EL100ED	Cat.: TL-022EL200ED	Cat.: TL-022EL400ED	Cat.: TL-022EL800ED
- TL-022EL01 Pancreatic Elastase Turbilatex [®] Reagent 1	1x24mL	2 x 24mL	4 x 24mL	8 x 24mL
- TL-022EL02 Pancreatic Elastase Turbilatex [®] Reagent 2	1 x 4mL	1 x 6mL	2 x 6mL	4 x 6mL
- TL-022EL70/ TL-022EL71/ TL-022EL72/ TL-022EL73/ TL-022EL74/ TL-022EL75 Pancreatic Elastase Turbilatex [®] Calibrators			6 x 1mL	
- TL-022EL08/ TL-022EL09 Pancreatic Elastase Turbilatex [®] Control 1&2		4 x 1mL		

Materials and instruments that could be required but no provided

- Specimen collection container.
- Disposable gloves and laboratory equipment.
- Automated analyser.
- MST-0019U Universal Turbilatex[®] Sample Collection Vial.

SPECIMEN COLLECTION AND PREPARATION

Collect enough quantity of human solid stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. If not immediately tested, freeze the stored samples at -20 °C maximum 6 months. In this case, the sample will be totally thawed and brought to room temperature (15-30°C) before testing. Homogenize stool samples as thoroughly as possible prior to preparation.

To process the collected stool samples using catalogue reference MST-0019U follow the indication of the instruction for use of the product.

ANALYTICAL PROCEDURE

R1 and R2 are ready to use.

Preparation of the calibration curve


For calibration use only the following materials:

TL-022EL70, TL-022EL71, TL-022EL72, TL-022EL73, TL-022EL74 and TL-022EL75, Pancreatic Elastase Turbilatex® Calibrator vials: liquid calibrator containing human pancreatic elastase at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a week is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

Pancreatic Elastase Turbilatex® Calibrator vials (calibrator 0 to 5) are ready to use.

Picture 1. Prepared calibrators for performing the calibration curve (5 points and blank)

Quality control


For quality control only use the following materials:

TL-022EL08 Pancreatic Elastase Turbilatex® Control 1 and TL-022EL09 Pancreatic Elastase Turbilatex® Control 2: liquid controls at two different concentrations of recombinant Pancreatic Elastase. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.

Pancreatic Elastase Turbilatex® Combo Control vials (Control 1 and Control 2) are ready to use.

Picture 2. Prepared controls C1 and C2 (check control concentration on the label)

INTERPRETATION OF RESULTS

Positive results: higher or equal than the cut-off fixed by the clinical lab.

Recommended: 100-200 µg of Pancreatic elastase E1/g of stool for diagnostic procedures.

Positive results determine the presence of normal level of Pancreatic Elastase in faeces.

Pancreatic Elastase Turbilatex® Combo cut-off value:

Pancreatic Elastase E1 values equal or higher than 200 µg of Pancreatic elastase E1/g of stool are indicative of a normal pancreatic exocrine functioning.

Pancreatic Elastase E1 values between 100 and 200 µg of Pancreatic elastase E1/g of stool are indicative of mild to moderate pancreatic exocrine insufficiency.

Pancreatic Elastase E1 values lower than 100 µg of Pancreatic elastase E1/g of stool are indicative of a severe pancreatic exocrine insufficiency.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Limit of detection (LOD): 1.07 µg hEL/g of stool. The lower limit of detection of **Pancreatic Elastase Turbilatex® Combo** was determined on 20 samples and 2 sample replicates as the mean value +2· SD.

Limit of quantification (LOQ): 5.3 µg hEL/g of stool. The lower limit of quantification is defined as the lowest actual amount of analysis that can be reliably detected; imprecision is < 20% as CV%.

Measuring range

Pancreatic Elastase Turbilatex® Combo assay measuring range is 5.3-625 µg hEL/g of stool. Samples more concentrated than 400 µg hEL/g of stool must be diluted for proper quantification by the user, using additional sample buffer.

Prozone

Studies have been made up to a concentration of 10000 µg hEL/g of stool and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Samples with concentrations up to 625 µg hEL/g of stool can be measured without inhibitory prozone effect.

Clinical sensitivity and specificity

The following tables provide information about the evaluations **Pancreatic Elastase Turbilatex® Combo**, has undergone with cut-off 200 µg hEL/g of stool:

N°	Evaluated test	Sample type
1	Schebo® Pancreatic Elastase 1, Schebo® Biotech AG	Human stool samples
2	BÜHLMANN fPELA® turbo, Bühlmann Diagnostics	Human stool samples
3	Schebo® Pancreatic Elastase 1, Schebo® Biotech AG	Human stool samples

Table 1 Resume of Clinical evaluations. Test and sample type.

N°	TP	FP	FN	TN	N	Sens. (%) CI (%)	Spec. (%) CI (%)	PPV (%) CI (%)	NPV (%) CI (%)	LR+ CI	LR- CI
1	10	2	0	28	40	100.0 (69.2-100.0)	93.3 (77.9-99.2)	83.3 (51.6-97.9)	100 (87.7-100.0)	11.836 (3.58-39.12)	0.049 (0.003-0.743)
2	22	0	3	56	81	88.0 (68.8-97.5)	100.0 (93.6-100.0)	100.0 (84.6-100.0)	94.9 (85.9-98.9)	96.65 (6.220-1565)	0.136 (0.051-0.360)
3	17	0	3	40	60	85.0 (62.1-96.8)	100.0 (91.2-100.0)	100.0 (80.5-100.0)	93.0 (80.9-98.5)	68.33 (4.318-1081)	0.169 (0.065-0.439)
Total	49	2	6	124	181	89.1 (77.8-95.9)	98.4 (94.4-99.8)	96.1 (86.5-99.5)	95.4 (90.2-98.3)	56.13 (14.15-222.7)	0.111 (0.052-0.236)

Table 2 Clinical evaluations and global results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI), Likelihood positive ratio (LR+), Likelihood negative ratio (LR-) for **Pancreatic Elastase Turbilatex® Combo**.

The results showed a high sensitivity and specificity to detect human pancreatic elastase (hEL) using **Pancreatic Elastase Turbilatex® Combo**.

Interferences

An evaluation was performed to determine the interferences of **Pancreatic Elastase Turbilatex® Combo**. No interferences were founded against other substances tested: BSA, Hemin, Ascorbic acid.

Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against porcine pancreatic elastase.

Precision

An evaluation was performed to determine the accuracy of **Pancreatic Elastase Turbilatex® Combo**. **Pancreatic Elastase Turbilatex® Combo** shows high in its reproducibility results on known concentration of human pancreatic elastase.

Linearity

Pancreatic Elastase Turbilatex® Combo using calibrator kit is linear in the calibration range of 18.8-400 µg hEL/g of stool.

The analytical performance of **Pancreatic Elastase Turbilatex® Combo** has been validated using the Tokyo Boeki Biolis 24i analyser. For more information about other compatible equipment, please refer to the corresponding application. The different turbidimetric equipment with which our turbidimetric products are compatible can be found in **Annex I**.

LIMITATIONS

1. **Pancreatic Elastase Turbilatex® Combo** should **only be used in human stool samples**. Proper faecal specimens must be obtained.
2. Values in the range between 100-200 µg/g should be consider as mild to moderate pancreatic insufficiency, and they should be reviewed by the specialist.
3. Conditions such as Crohn's disease, ulcerative colitis, and untreated celiac disease may cause elevated elastase levels regardless of pancreatic function, potentially resulting in falsely elevated readings (7).
4. Watery or loose stools may dilute elastase concentrations, leading to falsely low levels that could suggest pancreatic insufficiency (8).
5. In elderly patients and those with progressive pancreatic damage from cystic fibrosis, elastase levels may gradually decline, complicating the differentiation between progressive EPI and other conditions (1,2).
6. Long-term antibiotic use and other medications can alter gut microbiota, potentially affecting enzyme stability or elastase concentration in feces (7).
7. Sample characteristics like the presence of blood or high levels of fat can interfere with certain methods of elastase detection, potentially causing errors in interpretation (9).

EXPECTED VALUES

The prevalence of Exocrine pancreatic insufficiency (EPI) in the general population is unknown. It is most commonly associated with diseases of the exocrine pancreas, being a common manifestation in late stage chronic pancreatitis. This insufficiency can be measured indirectly by assessing quantitative changes in pancreatic secretion, e.g. of the proteolytic enzyme pancreatic elastase. (10).

Values of pancreatic elastase over 200 µg/g determine normal activity of the pancreas.

Values of pancreatic elastase between 100-200 µg/g determine mild to moderate pancreatic insufficiency and requires follow up.

Values of pancreatic elastase below 100 µg/g determine severe pancreatic insufficiency and requires additional diagnosis to implement the correct treatment (11).

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
1. Daftary A, Acton J, Heubi J, Amin R. **Fecal elastase-1: Utility in pancreatic function in cystic fibrosis.** Vol. 5, Journal of Cystic Fibrosis. 2006. p. 71–6.
2. Borowitz D, Baker SS, Duffy L, Baker RD, Fitzpatrick L, Gyamfi J, et al. **Use of fecal elastase-1 to classify pancreatic status in patients with cystic fibrosis.** Journal of Pediatrics. 2004 Sep;145(3):322–6.
3. Meyts I, Wuyts W, Proesmans M, De Boeck K. **Variability of fecal pancreatic elastase measurements in cystic fibrosis patients.** Vol. 1, Journal of Cystic Fibrosis. 2002.
4. Lindkvist B. **Diagnosis and treatment of pancreatic exocrine insufficiency.** World J Gastroenterol. 2013;19(42):7258–66.
5. Drucy Borowitz Ñ, Lin R, Baker SS. **Comparison of Monoclonal and Polyclonal ELISAs for Fecal Elastase in Patients With Cystic Fibrosis and Pancreatic Insufficiency** [Internet]. 2006. Available from: <http://journals.lww.com/jpgn>
6. Singh VK, Schwarzenberg SJ. **Pancreatic insufficiency in Cystic Fibrosis.** Journal of Cystic Fibrosis. 2017 Nov 1;16:S70–8.
7. Walkowiak J, Nousia-Arvanitakis S, Cade A, Kashirskaya N, Piotrowski R, Strzykala K, et al. **Fecal elastase-1 cut-off levels in the assessment of exocrine pancreatic function in cystic fibrosis.** Journal of Cystic Fibrosis. 2002 Dec 1;1(4):260–4.
8. Vanga RR, Tansel A, Sidiq S, El-Serag HB, Othman MO. **Diagnostic Performance of Measurement of Fecal Elastase-1 in Detection of Exocrine Pancreatic Insufficiency: Systematic Review and Meta-analysis.** Clin Gastroenterol Hepatol. 2018 Aug 1;16(8):1220-1228.e4.
9. Löser C, Möllgaard A, Folsch UR, Loser Mollgaard CA, Chr LUser PD. **Faecal elastase 1: a novel, highly sensitive, and specific tubeless pancreatic function test.** Vol. 39, Gut. 1996.
10. Capurso G, Traini M, Piciocchi M, Signoretti M, Arcidiacono PG. **Exocrine pancreatic insufficiency: prevalence, diagnosis, and management.** Clin Exp Gastroenterol. 2019;12:129.

11. Whitcomb DC, Buchner AM, Forsmark CE. **AGA Clinical Practice Update on the Epidemiology, Evaluation, and Management of Exocrine Pancreatic Insufficiency**: Expert Review. *Gastroenterology*. 2023;165:1292–301.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website: www.certest.es

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