# FOB TURBILATEX<sup>®</sup> Combo

## FOB latex turbidimetric assay

Cat. TL-022FB100ED, TL-022FB200ED, TL-022FB400ED, TL-022FB800ED

IU-TL - 022FB100ED-200ED-400ED-800ED en-es v.05

## CONTENT/CONTENIDO

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## ENGLISH

## INTENDED PURPOSE

**FOB Turbilatex® Combo** is a latex turbidimetric assay **for the quantitative detection of haemoglobin (faecal occult blood) in human stool samples.**

This assay is simple and widely applicable. This product is optimized for several automated analyser.

For professional *in vitro* diagnostic use only.

## SYNTHESIS

Colorectal cancer is the second leading cause of illness and death in Western world (1). The screening with faecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps (2). Colorectal cancer is also associated with local acute inflammatory reaction being visualized, in some cases, by white cell neutrophil scanning.

Haemoglobin is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases.

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site and amount of bleeding and the transit time through the gut (3). **FOB Turbilatex® Combo** test detect intact or nearly intact human haemoglobin, being a very specific technique for detecting loss of blood from the lower intestine.

## PRINCIPLE

**FOB Turbilatex® Combo** is a latex turbidimetric assay **for the quantitative detection of haemoglobin in human stool samples.** This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use of the test is to make a presumptive diagnosis of faecal occult blood (gastrointestinal bleeding).

FOB latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample it will react with the antibodies and form an aggregate. If no antigen is present in the

sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at determinate wavelenght and is proportional to the quantity of antigen contained in the sample.

## PRECAUTIONS

---

- For professional *in vitro* diagnostic use only.
- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in Turbidimetric technique, autoanalyzer use and *in vitro* diagnostic procedures.
- The analyser should be ready before performing any assay.
- Read the instructions for use carefully before use the kit. Do not perform the test without fully understanding the information on procedures, safely precautions and limitations described in this instruction for use.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Do not freeze to maintain the integrity and performance of the test.
- Do not mix screw caps from different reagent vials.
- The kit does not contain materials of human origin.
- Components provided in the kit are approved for use in **FOB Turbilatex® Combo**. Do not mix or use with components from other product kit or lots. Do not use with other commercial kit components (from another supplier). The test's satisfactory performance is ensured only when components from the same lot of **FOB Turbilatex® Combo** are used.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples.
- **FOB Turbilatex® Combo** do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their

declaration. Material Safety Data Sheet are not included with this device, nevertheless it could be consulted upon request to Certest Biotec.

- **FOB Turbilatex® Combo** has been validated with the analyser Biolis from Tokyo Boeki. For different analysers, Haemoglobin should be measured according to the specific application parameters for each specific chemistry analyzer. More information about the application notes available can be found on our website.
- Homogenise the reagents before use, make sure that the reagents are free of bubbles.
- Evaporation of calibrators and controls can lead to incorrect results. Follow the storage recommendations detailed in the following section.
- If the values obtained for the controls are not valid, repeat the quality control measurement with new controls. If the control values are still invalid, recalibrate the equipment.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- The solutions should be discarded in a proper container after testing following local regulations.

## STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). **Do not freeze.**

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight and conserved at 2-8°C.

## REAGENTS

### Materials provided

<u>References</u>	<u>Cat.:</u> <u>TL-022FB100ED</u>	<u>Cat.:</u> <u>TL-022FB200ED</u>	<u>Cat.:</u> <u>TL-022FB400ED</u>	<u>Cat.:</u> <u>TL-022FB800ED</u>
- TL-022FB01 FOB Turbilatex® Reagent 1	1x22mL	2 x 22mL	4 x 22mL	8 x 22mL
- TL-022FB02 FOB Turbilatex® Reagent 2	1 x 8mL	1 x 13mL	2 x 13mL	4 x 13mL
- TL-022FB70/ TL-022FB71/ TL- 022FB72- TL-022FB73/ TL- 022FB74/ TL-022FB75 FOB Turbilatex® Calibrators		6 x 1mL		
- TL-022FB08/ TL-022FB09 FOB Turbilatex® Control 1&2		4 x 1mL		

### Materials and instruments that could be required but no provided

- Specimen collection container.

- Disposable gloves and laboratory equipment.
- Automated analyser.
- MST-0019U Universal Turbilatex® Sample Collection Vial.

## SPECIMEN COLLECTION AND PREPARATION

Collect enough quantity of human stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 3 days prior to testing. If not immediately tested, freeze the stored samples at -20 °C maximum 6 months. In this case, the sample will be totally thawed and brought to room temperature (15-30°C) before testing. Homogenize stool samples as thoroughly as possible prior to preparation.

**To process the collected stool samples using catalogue reference MST-0019U follow the indication of the instruction for use of the product.**

## ANALYTICAL PROCEDURE

R1 and R2 are ready to use.

### Preparation of the calibration curve


For calibration use only the following materials:

TL-022FB70, TL-022FB71, TL-022FB72, TL-022FB73, TL-022FB74 and TL-022FB75, FOB Turbilatex® Calibrator vials: liquid calibrator containing human haemoglobin at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a week is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

FOB Turbilatex® Calibrator vials (calibrator 0 to 5) are ready to use.

*Picture 1.* Prepared calibrators for performing the calibration curve (5 points and blank)

## Quality control


For quality control only use the following materials:

TL-022FB08 FOB Turbilatex® Control 1 and TL-022FB09 FOB Turbilatex® Control 2: liquid controls at two different concentrations of human haemoglobin. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.

**FOB Turbilatex® Combo** Control vials (Control 1 and Control 2) are ready to use.

**Picture 2.** Prepared controls C1 and C2 (check control concentration on the label)

## INTERPRETATION OF RESULTS

It is recommended that each laboratory/end-user establish their own cut-off value based on specific clinical and demographic needs, as well as the operational requirements of their screening programs. Adjusting the cut-off enables programs to optimize the trade-off between detecting significant neoplasia and minimizing unnecessary colonoscopies.

Therefore, the following reference values are to be considered as recommendations.

Recommended cut-off value: 10 µg hHb/g in stool samples (4).

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity

The Limit of Detection (LOD) and Quantification (LOQ) were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A requirements.

Limit of detection (LOD): 32.33 ng hHb/mL.

Limit of quantification (LOQ): 34.60 ng hHb/mL. The lower limit of quantification is defined as the lowest actual amount of analysis that can be reliably detected. The Upper Limit of Quantification, it has not been determined since the LOQ experiment has proved that there is good quantification up to the 1000 ng hHb/mL point, which is the maximum point of the **FOB Turbilatex® Combo** calibration curve. This point can be quantified with a coefficient of variation lower than the %CV goal (13%).

## Measuring range

**FOB Turbilatex® Combo** assay measuring range is 34.60-1953 ng hHb/mL. Samples more concentrated than 1000 ng hHb/mL must be diluted for proper quantification by the user, using additional sample buffer.

## Prozone

Studies have been made up to a concentration of 1 mg of hHb/mL of stool and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Samples with concentrations up to 1953 ng hHb/mL can be measured without inhibitory prozone effect.

## Clinical sensitivity and specificity

An evaluation was performed comparing **FOB Turbilatex® Combo** against another commercially available turbidimetric assay, which was considered as a gold standard. This evaluation has undergone with cut-off 10 µg of hHb/g of stool:

	Evaluation criteria			Total
		+	-	
FOB Turbilatex® Combo	+	28	1	29
	-	4	46	50
	Total	32	47	79

Table 1 **FOB Turbilatex® Combo** performance versus Evaluation criteria.

FOB Turbilatex® vs Evaluation criteria		
	Mean Value	95% confidence interval
Sensitivity	87.5%	71.0-96.5%
Specificity	97.9%	88.7-99.9%
PPV	96.6%	82.2-99.9%
NPV	92.0%	80.8-97.8%
LR+	41.13	5.89-287.2

FOB Turbilatex® vs Evaluation criteria		
	Mean Value	95% confidence interval
LR-	0.128	0.051-0.320

Table 2 **FOB Turbilatex® Combo** clinical results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI), Likelihood positive ratio (LR+), Likelihood negative ratio (LR-) for **FOB Turbilatex® Combo**.

The results showed a high sensitivity and specificity to detect human haemoglobin (hHb) using **FOB Turbilatex® Combo**.

## Interferences

An evaluation was carried out to determine the interferences of **FOB Turbilatex® Combo** taking into account the recommendations of the CLSI EP07 and EP37 guidelines. No interferences were founded against the exogenous or endogenous substances tested:

Substance	Concentration
Acenocoumarol	1.20·10 <sup>-1</sup> mg/dL
Amitriptilin	4.80·10 <sup>-2</sup> mg/dL
Biotin	51·10 <sup>-1</sup> mg/dL
Ciprofloxacin	1.20 mg/dL
Diazepam	3.00 mg/dL
Macrogol 3350	3.00 mg/dL
Naproxen	3.60·10 <sup>1</sup> mg/dL
Omeoprazol (gastric protector)	8.40·10 <sup>-1</sup> mg/dL
Rifaximin	4.05·10 <sup>-4</sup> mg/dL
Oral saline solution	1.30·10 <sup>2</sup> mg/dL
Sulfasalazine	7.50 mg/dL
Barium sulfate	6.30·10 <sup>1</sup> mg/dL
Tetracycline Hydrochloride	1.20·10 <sup>-1</sup> mg/dL
Vancomycin	1.20·10 <sup>1</sup> mg/dL
Ascorbic acid	5.25 mg/dL
Glucose	1.00·10 <sup>3</sup> mg/dL
Triglycerides	1.00·10 <sup>2</sup> mg/dL
Hemin	1.00·10 <sup>-2</sup> mmol/L
Albumin	3.50·10 <sup>1</sup> g/L
Bilirubin	4.00·10 <sup>1</sup> mg/dL
Leukocytes	2.00·10 <sup>7</sup> cell/mL
Mucin	5.00%
Urea	4.29·10 <sup>1</sup> mmol/L
Uric acid	1.40·10 <sup>3</sup> µmol/L

Table 3 Possible interference substances evaluated.

## Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against other faecal markers occasionally present in faeces, such as:

human calprotectin, bovine and pig haemoglobin, human lactoferrin, bovine and human transferrin.

The cross reactions with bovine, pig haemoglobins and bovine transferrin are not present at concentrations below 0.20 µg/mL, 0.90 µg/mL and 250 µg/mL respectively. Take note that these concentrations are higher than the concentrations that could be found on the faecal samples.

## Precision

An evaluation was performed to determine the accuracy of **FOB Turbilatex® Combo**. Within-laboratory and repeatability were determined according to CLSI EP05 using a standardised study design of 80 replicates per sample were evaluated (5 days x 4 runs x 4 replicates) and with an acceptance criterion of 20% CV.

Sample	N	Mean (ng/mL)	Repeatability		Total	
			Sd	CV%	Sd	CV%
1	80	50.37	4.83	9.6%	5.62	11.15%
2	80	124.40	13.57	10.8%	18.62	14.87%
3	80	232.78	3.46	1.5%	11.46	4.92%
4	80	320.80	3.97	1.2%	18.44	5.75%
5	80	403.55	6.57	1.6%	43.52	10.78%
6	80	561.96	8.78	1.6%	13.16	2.34%
7	80	992.05	50.49	5.1%	51.22	5.16%

Table 4 Repeatability results for **FOB Turbilatex® Combo** including mean, standard deviation (Sd) and the coefficient of variation (CV%).

Similarly, following the CLSI EP05 recommendations for reproducibility analysis, 80 replicates per sample were evaluated with three lots (5 days x 4 runs x 4 replicates) and with an acceptance criterion of 20% CV.

Sample	N	Mean (ng/mL)	Repeatability		Within lot		Reproducibility	
			Sd	CV%	Sd	CV%	Sd	CV%
1	80	49.466	4.067	8.2%	7.111	14.4%	7.111	14.4%
2	80	124.443	8.666	7.0%	16.788	13.5%	17.241	13.9%
3	80	238.002	9.261	3.9%	15.706	6.6%	15.868	6.7%
4	80	323.423	2.936	0.9%	17.057	5.3%	17.372	5.4%
5	80	411.280	4.693	1.1%	39.011	9.5%	39.180	9.5%

Sample	N	Mean (ng/mL)	Repeatability		Within lot		Reproducibility	
			Sd	CV%	Sd	CV%	Sd	CV%
6	80	590.760	6.462	1.1%	29.908	5.1%	48.175	8.2%
7	80	1011.871	33.099	3.3%	43.607	4.3%	45.049	4.5%

Table 5 Reproducibility results for **FOB Turbilatex® Combo** including mean, standard deviation (Sd) and the coefficient of variation (CV%).

**FOB Turbilatex® Combo** shows high in its reproducibility results on known concentration of haemoglobin.

## Linearity

The linearity study was performed following CLSI EP06 using samples with assigned values. The allowable deviation from linearity (ADL) is  $\pm 10\%$  considering indications of EP06 CLSI guideline regarding ADL.

**FOB Turbilatex® Combo** using calibrator kit is linear in the calibration range of 40-1000 ng hHb/mL.

The analytical performance of **FOB Turbilatex® Combo** has been validated using the Tokyo Boeki Biolis 24i analyser. For more information about other compatible equipment, please refer to the corresponding application. The different turbidimetric equipment with which our turbidimetric products are compatible can be found in **Annex I**.

## LIMITATIONS

1. **FOB Turbilatex® Combo** should **only be used in human stool samples**. Proper faecal specimens must be obtained.
2. Positive results determine the presence of human haemoglobin in faecal samples. A positive result should be followed up with additional diagnostic invasive procedures to determine the exact cause and source of the blood in the stool.
3. If symptoms or situation still persist, haemoglobin determination should be carried out invasive techniques. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.
4. This test may be less sensitive for detecting upper gastrointestinal bleeding because blood degrades as it passes through the gastrointestinal tract.
5. Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.

## EXPECTED VALUES

---

Colorectal cancer (CRC) is a leading cause of cancer incidence and mortality globally. According to the World Health Organization, it is the third most common cancer worldwide, accounting for approximately 10% of all cancer cases and is the second leading cause of cancer-related deaths worldwide (1). In Europe, it ranks as the second most common cancer after breast cancer, with over 500,000 new cases reported annually (5). The increasing incidence among individuals under 50 has prompted a shift in screening recommendations, with guidelines now advocating routine screening starting at age 45 (6).

Fecal immunochemical tests (FIT) are commonly employed as a non-invasive screening method, detecting hemoglobin (hHb) in stool samples. Globally, enhancing early detection through these screening strategies is crucial, as early-stage CRC is more treatable, leading to better patient outcomes (2). Implementing standardized screening protocols can significantly reduce CRC-related morbidity and mortality across diverse populations.

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

1. World Health Organization. **Colorectal cancer. Fact sheet.** [Internet]. 2023 [cited 2024 Oct 29]. Available from: <https://www.who.int/news-room/fact-sheets/detail/colorectal-cancer>
2. Young GP, Symonds EL, Allison JE, Cole SR, Fraser CG, Halloran SP, et al. **Advances in Fecal Occult Blood Tests: The FIT Revolution.** Dig Dis Sci. 2015 Mar 1;60(3):609–22.
3. Ransohoff DF, Lang CA. **Screening for Colorectal Cancer with the Fecal Occult Blood Test: A Background Paper.** 1997;
4. NICE. **Quantitative faecal immunochemical testing to guide colorectal cancer pathway referral in primary care** [Internet]. 2023 Aug [cited 2024 Nov 26]. Available from: [www.nice.org.uk/guidance/dg56](http://www.nice.org.uk/guidance/dg56)
5. International Agency for Research on Cancer. **Cancer Today. Absolute numbers, Incidence, Both sexes, in 2022** [Internet]. 2022 [cited 2024 Oct 29]. Available from: [https://gco.iarc.who.int/today/en/dataviz/pie?mode=cancer&group\\_populations=1&populations=908](https://gco.iarc.who.int/today/en/dataviz/pie?mode=cancer&group_populations=1&populations=908)
6. American Cancer Society. **Colorectal Cancer Early Detection, Diagnosis, and Staging** [Internet]. 2024. Available from: [https://seer.cancer.gov/csr/1975\\_2016/](https://seer.cancer.gov/csr/1975_2016/),

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website: [www.certest.es](http://www.certest.es)

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