CerTest Rapid Test

CerTest
BIOTEC

CerTest Strep A

ONE STEP *Group A Streptococcus* CARD TEST

CE IVD

The test should be used by professional in IVD devices

These instructions for use apply to the following references:

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T820201PC; T820201PCC

Initial part of the Instructions for use is common for all the references/

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Specific parts for each of the references are described in the corresponding annexes (see tables below)/

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KIT 20 CARD GROUP A STREPTOCOCCUS WITH POSITIVE CONTROL (SEE ANNEX 1) /

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KIT 20 CARD GROUP A STREPTOCOCCUS WITH POSITIVE AND NEGATIVE CONTROL (SEE ANNEX 2) /

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PRODUCT / PRODUCTO	REFERENCE / REFERENCIAS
CerTest Strep A Card (20 tests kit with swabs + PC + NC)	T820201PCC

Table A 2. Reference for Kit 20 card of CerTest Strep A with positive and negative control / Referencia para kit de 20 carcassas de CerTest Strep A con controles positivo y negativo.

1. Quick Guide/Guía rápida

Quick Guide

CerTest Strep A

CerTest
BIOTEC

1

Ref. T820201PC

Ref. T820201PCC

C+

C-

88:88:88

(x20) (x20) (x1) (x1) (x20)

2

1°

2°

(x5) (x5)

(1 min)

3

x3

4

88:10:00

10'

C T

C T

- +

5

Content/ Contenido

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1.	Quick Guide/Guía rápida	3
----	-------------------------------	---

ENGLISH

2.	Intended use	7
3.	Summary and Explanation	7
4.	Test principle	7
5.	Storage and stability	8
6.	Precautions for users	8
7.	Specimen collection and preparation	10
8.	Limitations of the test	11
9.	Expected values	12
10.	Performance characteristics	12
A.	Analytical sensitivity (detection limit)	12
B.	Clinical sensitivity and specificity	12
C.	Hook Effect	13
D.	Cross Reactivity	13
E.	Interferences	14
F.	Repeatability and Reproducibility	15
	ANNEX 1	16
I.	Materials	16
II.	Test procedures (samples)	16
III.	Interpretation of results	17
IV.	Quality Control	18
	A. External control. CerTest Strep A Positive Control.....	18
A.1	Positive control procedure	18
A.1.1	Interpretation positive control procedure	20
	ANNEX 2	21
I.	Materials	21
II.	Test procedures (samples)	21
III.	Interpretation of results	22
IV.	Quality Control	23
	A. External control. CerTest Strep A Positive Control.....	23

A.1	Positive control procedure.....	23
A.1.1	Interpretation positive control procedure	25
B.	External negative control. CerTest Negative Control.....	25
B.1	Negative control procedure	25
B.1.1.	Interpretation negative control procedure	27

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ENGLISH

2. Intended use

CerTest Strep A one step card test is a coloured chromatographic immunoassay for the qualitative detection of Group A *Streptococcal* from throat swabs and suspected Group A *Streptococcal* colonies recovered from culture.

CerTest Strep A card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of Group A *Streptococcal* respiratory infection.

3. Summary and Explanation

Group A *Streptococcal* infections are caused by Group A *Streptococcus* (GAS), a gram-positive bacterium responsible for a variety of health problems (1). These infections can range from mild human infections such as pharyngitis and impetigo to severe life-threatening conditions such as necrotizing fasciitis, often referred as "flesh-eating disease", myositis or streptococcal toxic shock syndrome (STSS) (1,2).

GAS infection most frequently occurs in children between the ages of 5 and 15 years and primarily in the cooler months of winter and early spring. Patients typically present with a sudden-onset sore throat, pain upon swallowing and fever. Without treatment, symptoms may last 8 to 10 days and patients are infectious for up to 1 week following the acute presentation (2,3).

Bacterial culture is considered the gold standard for the diagnosis of GAS due to high sensitivity and specificity. However, the incubation, isolation and subsequently identification of viable pathogens by this technique requires 24 to 48 hour or longer. Compared with throat swabs culture, rapid antigen detection tests (RADTs) offer diagnosis 10 minutes (4).

For the accurate diagnosis of GAS, the Infectious Disease Society of America (IDSA) recommends using a rapid test in patients with a modest probability of GAS infection, treating those with a positive rapid test and withholding antibiotics in rapid test negative patients. Negative results should be backed up with bacterial culture in pediatric patients and treating patients having positive cultures (3).

Rapid antigen testing to detect Group A *Streptococcal* (GAS) infection provides important information for the antibiotic decision-making patients presenting with acute pharyngitis and reduces the severity of symptoms and further complications (2).

4. Test principle

CerTest Strep A is based on the principle of a qualitative immunochromatographic assay for determination of Group A *Streptococcus* in human throat samples from patients suspected of Group A *Streptococcal* infection.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Group A *Streptococcus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti- Group A *Streptococcus*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to blue polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigen of the diluted sample reacts with the red-coloured conjugate complex (anti-Strep A monoclonal antibodies-red polystyrene microspheres) which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-Strep A antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no Group A *Streptococcus* antigens presence and yet, the antigen may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-Strep A antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control blue-conjugate complex and the control line will always appear. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

5. Storage and stability

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Once the **Reagent A** and **B** are opened, they must be kept in the conditions indicated on its label. Under these storage conditions, the reagents maintain their performance until expiry date indicated on its label. Do not freeze. Before performance the test, allow all the components to reach room temperature (15-30°C) prior to testing.

6. Precautions for users

- **For professional *in vitro* diagnostic use only.**
- Read the instructions for use carefully before use the kit.
- Do not use after expiration date.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use product if the protective external box or the protective aluminium pouches are opened or damaged upon arrival.

- Do not use the **Reagent A** and the **Reagent B** if they are opened or damaged upon arrival.
- Do not interchange caps between **Reagent A** and **Reagent B**.
- Do not use the test (and positive/negative controls if the final reference contains them) if desiccant material is not present or broken inside the aluminium pouch.
- Components provided in the kit are approved for use in **CerTest Strep A**. Do not mix or use with components from other product kit or lots. Do not use with other commercial kit components (from another supplier).
- Excepting from the **Reagent A** and the **Reagent B**, all the components of the kit are for single use. Do not reuse.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves or other personal protective equipment such as googles and mask that will be considered necessary. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous, as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples. Each sample must be correctly and unequivocally identified, in order to guarantee the correct traceability of the samples.
- A new test must be used for each sample to avoid contaminations errors.
- Clean up spills thoroughly using an appropriate disinfectant.
- Reagents contain preservatives (12-14% sodium nitrite). Avoid any contact with skin or mucous membrane. In accordance with Regulation (EC) No 1907/2006 (REACH), **CerTest Strep A** contains substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration. Material Safety Data Sheet are not included with this device, nevertheless it could be consulted upon request to CerTest Biotec.
- The certificate of analysis is not included with the device, however, it could be requested to Certest Biotec in case of need.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Visual interpretation of the results by coloured lines, interpretation of the results must be done by professional user without problems in colour interpretation.
- The tests (devices, testing tubes, disposable pipettes, positive/negative control if any) should be discarded in a proper biohazard container after testing. These containers should be discarded in accordance with local or national laws or regulations.

7. Specimen collection and preparation

Samples should be collected in a clean container and processed as soon as possible after collection. Follow proper infection control practices.


Samples can be preserved at room temperature if they are going to be analysed into the same day after taking them. The recommendation is processed the sample immediately after collection.

If this is not possible, throat samples could be stored in the refrigerator (2-8°C) for 8 hours prior for testing.

Specimen preparation (see illustration):

Throat swab method:


1. Remove the **swab** from its packaging.
2. Collect the specimen with the **swab** from the tonsils and/or the back of the throat. Be careful not to touch the teeth, tongue or cheek surfaces to avoid contamination with normal flora.
3. To ensure a good sample rotate the **swab** when sampling. Group A *Streptococcus* infects tonsilles, remember to swab on both sides.
4. Take the **swab** out of the mouth and process the throat swab specimen as soon as possible.

Suspected Group A *Streptococcal* colonies:

The Blood Agar is the method of choice for the isolation of Group A *Streptococcus* (anaerobic atmosphere, 72 hours/37°C). After 72 hours of incubation in blood the typical Group A *Streptococcal* colonies will be translucent white showing beta haemolysis.

1. Examine Blood plates after incubation period. Select *Streptococcus* typical colonies. Use a swab to pick up 1 or 3 colonies suspected (showing characteristic beta haemolysis).
2. Follow the instructions in the test procedure section to test the swab.

Continue immediately with the test procedure. See Annexes. Each sample vial (containing the sample) and cassette with which will be tested that sample, should be perfectly recognizable including patient's identification.

8. Limitations of the test

- 1) The test must be carried out within 15 minutes after opening the sealed bag.
- 2) The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- 3) This test does not differentiate between carriers and acute infection. Pharyngitis may be caused by organisms other than Group A *Streptococcus*.
- 4) **CerTest Strep A** should be used only with throat swabs or colonies taken directly from a plate. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper throat swab specimens must be obtained.
- 5) Positive results determine the presence of Group A *Streptococcus* in throat samples; nevertheless, a positive result should be followed up with additional laboratory techniques to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- 6) A negative result is not meaningful because of it is possible the antigens concentration in the throat samples is lower than the detection limit value. If the symptoms or situation persist, a Group A *Streptococcus* determination should be carried out, on a sample from an enrichment culture.
- 7) Excess blood or mucous on the swab specimen may interfere with test performance and may yield a false positive result. Avoid touching the tongue, cheeks and teeth and any bleeding areas of the mouth with the swab when collecting specimens.

9. Expected values

The World Health Organization (WHO) rank Group A Streptococcus (GAS) as the ninth leading infectious cause of human mortality (2).

GAS is a human-specific pathogen that is highly prevalent worldwide causing ~750 million infections per year, mostly at the throat and skin (5).

GAS accounts for 20% to 40% of cases of pharyngitis in children and 5% to 15% in adults. Most cases are benign and self-limiting within a week but suppurative complications or non-suppurative post-streptococcal diseases can occur (4).

10. Performance characteristics

A. Analytical sensitivity (detection limit)

Detection limit value (LoD) of **CerTest Strep A** is $9.8 \cdot 10^3$ CFU/mL of *Group A Streptococcus*. There is not any external standard approved as a general and reference standard. The internal standard used for validation of the CerTest products is inactivated Strep A antigen.

B. Clinical sensitivity and specificity

An external retrospective evaluation with throat specimens from patients with pharyngitis symptoms was performed using **CerTest Strep A** versus bacterial culture method (gold-standard). This is the evaluation criteria. The results were as follows:

		Evaluation criteria		
		+	-	Total
IC test: CerTest Strep A	+	97	4	101
	-	19	94	113
	Total	116	98	214

Table 1 CerTest Strep A performance versus evaluation criteria.

CerTest Strep A vs Evaluation criteria		
	Mean Value	95% confidence interval
Sensitivity	83.6%	75.6-89.8%
Specificity	95.9%	89.9-98.9%
PPV	96.0%	90.2-98.9%
NPV	83.2%	75.0-89.6%

Table 2 CerTest Strep A clinical results.

The most recent external retrospective evaluation performed to the test showed the following results: sensitivity 83.6% (95% confidence interval: 75.6-89.9%) and specificity: 95.9% (95% confidence interval: 89.9-98.9%).

The results showed a high sensitivity and specificity to detect *Group A Streptococcus* using **CerTest Strep A**.

The following tables provide information about the evaluations **CerTest Strep A** has undergone.

N°	Site	Sample type	Reference	Target
1	Al-Imanein Al-Kadhimein Medical City and Al-Zahra Primary Health Care Center (Bagdad, Iraq) November 2015-May 2016	Human throat samples	(17.1.1)	<i>Group A Streptococcus</i>
2	Public hospital (Madrid, Spain) 2008	Human throat samples	(17.1.2)	<i>Group A Streptococcus</i>
3	In-hospital NABL certified laboratory (India) December 2014- June 2015	Human throat samples	(17.1.3)	<i>Group A Streptococcus</i>
4	Public hospital (France) 2013	Human throat samples	(17.1.4)	<i>Group A Streptococcus</i>

Table 3 Resume of Clinical evaluations. Site, sample type, target.

N°	TP	FP	FN	TN	N	Sens. (%) CI (%)	Spec. (%) CI (%)	PPV (%) CI (%)	NPV (%) CI (%)
1	97	4	19	94	214	83.6 (75.6-89.8)	95.9 (89.9-98.9)	96.0 (90.2-98.9)	83.2 (75.0-89.6)
2	4	0	0	20	24	100.0 (39.8-100.0)	100.0 (83.2-100.0)	100.0 (39.8-100.0)	100.0 (83.2-100.0)
3	22	4	1	63	90	95.7 (78.1-100.0)	94.0 (85.4-98.3)	84.6 (65.1-95.6)	98.4 (91.6-100.0)
4	8	1	0	29	38	100.0 (63.1-100.0)	96.7 (82.8-99.9)	88.9 (51.8-99.7)	100.0 (88.1-100.0)
Total	131	9	20	206	366	86.8 (80.3-91.7)	95.8 (92.2-98.1)	93.6 (88.1-97.0)	91.2 (86.7-94.5)

Table 4 Clinical evaluations and global results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Positive Predictive Values (PPV), Negative Predictive Values (NPV), 95% interval of confidence (CI) for CerTest Strep A.

C. Hook Effect

CerTest Strep A does not show full inhibitory hook effect at concentration higher than $10^5 \times \text{LoD}$.

D. Cross Reactivity

An evaluation was performed to determine the cross reactivity of **CerTest Strep A** no cross reactivity against organism, pathogens, substances that could cause infections:

Cross reactivity										
Virus	Adenovirus	-	Coronavirus	-	Influenza B	-	Norovirus GII	-	Rotavirus	-
	Astrovirus	-	Influenza A	-	Norovirus GI	-	Respiratory Syncytial Virus	-		
Bacteria	Campylobacter coli	-	Clostridium perfringens	-	Legionella pneumophila	-	Salmonella typhimurium	-	Streptococcus pneumococcal	-
	Campylobacter jejuni	-	Escherichia coli O:157	-	Listeria monocytogenes	-	Shigella boydii	-	Yersinia enterocolitica O:3	-
	Clostridium difficile antigen GDH	-	Escherichia coli O:026	-	Salmonella enteritidis	-	Shigella dysenteriae	-	Yersinia enterocolitica O:9	-
	Clostridium difficile Toxin A	-	Escherichia coli O:111	-	Salmonella paratyphi	-	Shigella flexneri	-		
	Clostridium difficile Toxin B	-	Helicobacter pylori	-	Salmonella typhi	-	Shigella sonnei	-		
Parasites	Cryptosporidium Parvum	-	Entamoeba histolytica	-	Giardia lamblia	-				
Others	Pig haemoglobin	-	Bovine haemoglobin	-						

Table 5 Cross reactivity results.

E. Interferences

An evaluation was performed to determine the possible interferences of **CerTest Strep A**; no interferences against the substances tested were detected:

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metronidazole	3.0 mg/mL	No interference	Loperamide hydrochloride (Fortasec)	0.15 mg/mL	No interference
	Ampicillin	3.0 mg/mL	No interference	Heparin (Hibor)	350.0 IU/mL	No interference
	Oseltamivir	0.3 mg/mL	No interference	Almagato (Almax)	3.0 mg/mL	No interference
	Amantadine	0.3 mg/mL	No interference	Fosfamycin (Monurof)	3.0 · 10 ⁻³ mg/mL	No interference
	Ribavirin	3.0 mg/mL	No interference	Acetylcysteine (Fluimucil)	3.0 mg/mL	No interference
	Codeine (Toseina)	0.2 mg/mL	No interference	Dexketoprofen trometamol (Enantyum)	0.3 mg/mL	No interference
	Benzocaine (Angileptol)	0.3 mg/mL	No interference	Levofloxacin	3.0 mg/mL	No interference
	Cloperastine (Flutox)	0.3 mg/mL	No interference	Ciprofloxacin	0.3 mg/mL	No interference
	Carbocisteine (Iniston mucolítico)	3.0 mg/mL	No interference	Rifampicin (Rifaldin)	0.3 mg/mL	No interference

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
	Loratadine	0.3 mg/mL	No interference	Phenoxyethylpenicillin potassium	3.0 mg/mL	No interference
	Dexchlorpheniramine (Polaramine)	0.3 mg/mL	No interference	Ambroxol hydrochloride (Mucosan)	0.3 mg/mL	No interference
	Ebastine (Ebastel)	3.0 mg/mL	No interference	Macrogol 3350 (Movicol)	3.0 mg/mL	No interference
	Acetyl Salicylic (Adiro)	0.3 mg/mL	No interference	Lysine Carbocysteinate (Pectox)	3.0 · 10 ⁻² mg/mL	No interference
	Ibuprofen (Espidifen)	3.0 mg/mL	No interference	Hydroxyzine dihydrochloride	0.3 mg/mL	No interference
	Paracetamol (Dolocatil)	5.0 mg/mL	No interference	Lorazepam	3.0 · 10 ⁻³ mg/mL	No interference
	Metamizole (Nolofil)	5.0 mg/mL	No interference	Amoxicillin	3.0 mg/mL	No interference
	Prednisone	0.3 mg/mL	No interference	Mercaptopurine	0.3 mg/mL	No interference
	Omeprazole	2.0 · 10 ⁻³ mg/mL	No interference	Biotine	0.1 mg/mL	No interference
	Naso GEL	0.9 mg/mL	No interference	Sore Throat Phenol spray	0.5 mg/mL	No interference
	CVS Nasal Spray (Cromolyn)	4 mg/mL	No interference	Tobramycin	0.3 mg/mL	No interference
	Afrin (Oxymetazoline)	5.0 · 10 ⁻² mg/mL	No interference	Mupirocin	2.5 · 10 ⁻² mg/mL	No interference
	CVS Nasal Drops (Phenylephrine)	10.0 mg/mL	No interference	Fluticasone Propionate	5.0 · 10 ⁻² mg/mL	No interference
	ZICAM	1.0 mg/mL	No interference	Homeopathic	DIL 1/10	No interference
	Phenylpropanolamine	0.3 mg/mL	No interference			
Endogenous	Human haemoglobin	5.0 µg/mL	No interference	Human lactoferrin	5.0 µg/mL	No interference
	Human transferrin	0.5 µg/mL	No interference	Mucine	0.5 %	No interference
	Human calprotectin	5.0 µg/mL	No interference	Human blood	5.0 %	No interference

Table 6 Interferences results.

F. Repeatability and Reproducibility

Following the study of repeatability and reproducibility of **CerTest Strep A** performed using different internal samples, negative and positive, no differences have been observed within the evaluations.

ANNEX 1

KIT 20 CARD GROUP A STREPTOCOCCUS. Catalogue reference: T820201PC

Annex for the following reference:

PRODUCT	REFERENCE
CerTest Strep A Card (20 tests kit with swabs + PC)	T820201PC

Table A 1. Reference for kit 20 card of CerTest Strep A with positive control.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
-20x CerTest Strep A (card test format)	- Disposable gloves (other personal protective equipment that will be considered necessary)
-1x Reagent A	- Timer
-1x Reagent B	- Specimen collection container
-20x Swab	
-20x Testing tube	
-20x Disposable pipette	
-1x CerTest Strep A Positive Control	
- 1x IFU Instructions for use (including Quick guide use).	


Table A 1.1 Materials necessities for test performance.

II. Test procedures (samples)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to use a separate **disposable pipette** and **testing tube** and device for each sample.
2. Add 5 drops (1) of **Reagent A** (light pink) and 5 drops (2) of **Reagent B** in a **testing tube** (the solution should turn light yellow (colourless)) and immediately introduce the **swab** into the **testing tube**.
3. Mix the solution by rotating the **swab** forcefully against the side of the **testing tube** at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (3). Extract as much liquid as possible from the **swab**, squeezing the sides of the **testing tube** or rotating the **swab** against the side of the **testing tube** as the **swab** is withdrawn. Discard the **swab**.
4. Remove **CerTest Strep A** from its sealed bag just before using it. The test must be carried out immediately after opening the sealed bag. Do not performed the test more than 15 minutes after opening the bag.

5. Dispense 3 drops from the **testing tube**, in the circular window marked with the letter S (4).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

III. Interpretation of results

(please, refer to the illustration below)


	CerTest Strep A	Interpretation of the results
1.	NEGATIVE	 <p>There is no Group A <i>Streptococcus</i> (<i>Strep A</i>) presence in the sample. No infection caused by Group A <i>Streptococcus</i>. A negative result should be followed up with additional laboratory techniques to confirm the results.</p>
	BLUE	
2.	POSITIVE	 <p>There is Group A <i>Streptococcus</i> (<i>Strep A</i>) presence in the sample. Possible infection caused by Group A <i>Streptococcus</i>. The presence of these lines (control blue line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.</p>
	BLUE-RED	
3.	INVALID	 <p>Invalid result, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.</p>
	Any other result	

Table A 1.2 CerTest Strep A performance interpretation of the results.

NEGATIVE: Only one BLUE line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the BLUE control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (BLUE) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control

CerTest Strep A (card test) presents two controls:

INTERNAL CONTROL: control line

Internal control: internal procedural control is included in the test. A blue line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROL: POSITIVE CONTROL SWAB

External control: CerTest Strep A Positive Control (positive control swab).

A. External control. CerTest Strep A Positive Control

CerTest Strep A Positive Control is an external quality control for **CerTest Strep A**. Use **CerTest Strep A Positive Control** to check that the extraction reagents and the test are working properly.


Store as packaged in the sealed pouch at 2-30°C. **CerTest Strep A Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest Strep A Positive Control** must remain in the sealed pouch until use. Do not freeze.

A.1 Positive control procedure

Allow tests and CerTest Strep A Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Proceed to use a separate **disposable pipette** and **testing tube** and device for each sample.
2. Add 5 drops (1) of **Reagent A** (light pink) and 5 drops (2) of **Reagent B** in a **testing tube** (the solution should turn light yellow (colourless)) and immediately introduce **CerTest Strep A Positive Control** into the **testing tube**.

3. Mix the solution by rotating **CerTest Strep A Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Best results are obtained when the specimen is vigorously extracted in the solution (3). Extract as much liquid as possible from **CerTest Strep A Positive Control**, squeezing the sides of the **testing tube** or rotating **CerTest Strep A Positive Control** against the side of the **testing tube** as **CerTest Strep A Positive Control** is withdrawn. Discard **CerTest Strep A Positive Control**.
4. Remove **CerTest Strep A** from its sealed bag just before using it. The test must be carried out immediately after opening the sealed bag. Do not performed the test more than 15 minutes after opening the bag.
5. Dispense 3 drops from the **testing tube**, in the circular window marked with the letter S (4).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Positive controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1.1 Interpretation positive control procedure



CerTest Strep A Positive Control			Interpretation of the results
1.	POSITIVE		A BLUE line (control line (C)) and a RED line (test line (T)) appear across the results window during the test performance due to presence of Group A Streptococcus in CerTest Strep A Positive Control .
	BLUE-RED		
2.	INVALID		Total absence of the control coloured line (BLUE), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (BLUE). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new CerTest Strep A Positive Control . If the situation still persists, discontinue using the test kit and contact your local distributor. Any result read after more than 10 minutes will be considered invalid.
	Any other result		

Table A 1.3 CerTest Strep A Positive Control performance interpretation of the results.

ANNEX 2

KIT 20 CARD GROUP A STREPTOCOCCUS. Catalogue reference: T820201PCC

Annex for the following reference:

PRODUCT	REFERENCE
CerTest Strep A Card (20 tests kit with swabs + PC + NC)	T820201PCC

Table A 2. Reference for kit 20 card of CerTest Strep A with positive and negative control.

I. Materials

MATERIALS PROVIDED	MATERIALS REQUIRED BUT NO PROVIDED
-20x CerTest Strep A (card test format) -1x Reagent A -1x Reagent B -20x Swab -20x Testing tube -20x Disposable pipette -1x CerTest Strep A Positive Control -1x CerTest Negative Control - 1x IFU Instructions for use (including Quick guide use).	- Disposable gloves (other personal protective equipment that will be considered necessary) - Timer - Specimen collection container


Table A 2.1 Materials necessities for test performance.

II. Test procedures (samples)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to use a separate **disposable pipette** and **testing tube** and device for each sample.
2. Add 5 drops (1) of **Reagent A** (light pink) and 5 drops (2) of **Reagent B** in a **testing tube** (the solution should turn light yellow (colourless)) and immediately introduce the **swab** into the **testing tube**.
3. Mix the solution by rotating the **swab** forcefully against the side of the **testing tube** at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (3). Extract as much liquid as possible from the **swab**, squeezing the sides of the **testing tube** or rotating the **swab** against the side of the **testing tube** as the **swab** is withdrawn. Discard the **swab**.

4. Remove **CerTest Strep A** from its sealed bag just before using it. The test must be carried out immediately after opening the sealed bag. Do not performed the test more than 15 minutes after opening the bag.
5. Dispense 3 drops from the **testing tube**, in the circular window marked with the letter S (4).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

III. Interpretation of results

(please, refer to the illustration below)

CerTest Strep A		Interpretation of the results	
1.	NEGATIVE		There is no Group A <i>Streptococcus</i> (<i>Strep A</i>) presence in the sample. No infection caused by Group A <i>Streptococcus</i> . A negative result should be followed up with additional laboratory techniques to confirm the results.
	BLUE		
2.	POSITIVE		There is Group A <i>Streptococcus</i> (<i>Strep A</i>) presence in the sample. Possible infection caused by Group A <i>Streptococcus</i> . The presence of these lines (control blue line and test red line), independent of their intensity should be considered as a positive result. A positive result should be followed up with additional laboratory techniques to confirm the results.
	BLUE-RED		
3.	INVALID		Invalid result, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid.
	Any other result		

Table A.2.2 CerTest Strep A performance interpretation of the results.

NEGATIVE: Only one BLUE line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the BLUE control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (BLUE) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

IV. Quality Control

CerTest Strep A (card test) presents two controls:

INTERNAL CONTROL: control line

Internal control: internal procedural control is included in the test. A blue line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

EXTERNAL CONTROL: POSITIVE CONTROL SWAB

External control: CerTest Strep A Positive Control (positive control swab).

A. External control. CerTest Strep A Positive Control

CerTest Strep A Positive Control is an external quality control for **CerTest Strep A**. Use **CerTest Strep A Positive Control** to check that the extraction reagents and the test are working properly.


Store as packaged in the sealed pouch at 2-30°C. **CerTest Strep A Positive Control** is stable through the expiration date printed on the sealed pouch. **CerTest Strep A Positive Control** must remain in the sealed pouch until use. Do not freeze.

A.1 Positive control procedure

Allow tests and CerTest Strep A Positive Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Proceed to use a separate **disposable pipette** and **testing tube** and device for each sample.

2. Add 5 drops (1) of **Reagent A** (light pink) and 5 drops (2) of **Reagent B** in a **testing tube** (the solution should turn light yellow (colourless)) and immediately introduce **CerTest Strep A Positive Control** into the **testing tube**.
3. Mix the solution by rotating **CerTest Strep A Positive Control** forcefully against the side of the **testing tube** at least 15 seconds. Best results are obtained when the specimen is vigorously extracted in the solution (3). Extract as much liquid as possible from **CerTest Strep A Positive Control**, squeezing the sides of the **testing tube** or rotating **CerTest Strep A Positive Control** against the side of the **testing tube** as **CerTest Strep A Positive Control** is withdrawn. Discard **CerTest Strep A Positive Control**.
4. Remove **CerTest Strep A** from its sealed bag just before using it. The test must be carried out immediately after opening the sealed bag. Do not performed the test more than 15 minutes after opening the bag.
5. Dispense 3 drops from the **testing tube**, in the circular window marked with the letter S (4).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Positive controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

A.1.1 Interpretation positive control procedure



CerTest Strep A Positive Control			Interpretation of the results
1.	POSITIVE		<p>A BLUE line (control line (C)) and a RED line (test line (T)) appear across the results window during the test performance due to presence of Group A <i>Streptococcus</i> in CerTest Strep A Positive Control.</p>
	BLUE-RED		
2.	INVALID		<p>Total absence of the control coloured line (BLUE), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (BLUE). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new CerTest Strep A Positive Control.</p> <p>If the situation still persists, discontinue using the test kit and contact your local distributor.</p> <p>Any result read after more than 10 minutes will be considered invalid.</p>
	Any other result		

Table A 2.3 CerTest Strep A Positive Control performance interpretation of the results.

B. External negative control. CerTest Negative Control

CerTest Negative Control is an external quality control for **CerTest Strep A** (card test). Use **CerTest Negative Control** to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C. **CerTest Negative Control** is stable through the expiration date printed on the sealed pouch. **CerTest Negative Control** must remain in the sealed pouch until use. Do not freeze.


B.1 Negative control procedure

Allow tests and CerTest Negative Control to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Proceed to use a separate **disposable pipette** and **testing tube** and device for each sample.
2. Add 5 drops (1) of **Reagent A** (light pink) and 5 drops (2) of **Reagent B** in a **testing tube** (the solution should turn light yellow (colourless)) and immediately introduce **CerTest Negative Control** into the **testing tube**.
3. Mix the solution by rotating **CerTest Negative Control** forcefully against the side of the **testing tube** at least 15 seconds. Best results are obtained when the specimen is vigorously extracted in the solution (3). Extract as much liquid as possible from **CerTest Negative Control**, squeezing the sides of the **testing tube**

or rotating **CerTest Negative Control** against the side of the **testing tube** as **CerTest Negative Control** is withdrawn. Discard **CerTest Negative Control**.

4. Remove **CerTest Strep A** from its sealed bag just before using it. The test must be carried out immediately after opening the sealed bag. Do not performed the test more than 15 minutes after opening the bag.
5. Dispense 3 drops from the **testing tube**, in the circular window marked with the letter S (4).
6. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Negative controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

B.1.1. Interpretation negative control procedure

CerTest Negative Control		Interpretation of the results
1. VALID RESULT		Only a BLUE line (control line (C)) appear across the results window during the test performance.
2. INVALID		<p>Total absence of the control coloured line (BLUE), regardless the appearance or not of the test line (RED), or present of both lines (RED test and BLUE control lines). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new CerTest Negative Control.</p> <p>If the situation still persists, discontinue using the test kit and contact your local distributor.</p> <p>Any result read after more than 10 minutes will be considered invalid</p>

Table A 2.4 CerTest Negative Control performance interpretation of the results.

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11. Reference/Bibliografía

1. Weiss KA LM. **Group A Streptococcus invasive infections: a review.** *Can J Chem.* 1997;40(1):18–25.
2. Walker MJ, Barnett TC, McArthur JD, Cole JN, Gillen CM, Henningham A, et al. **Disease manifestations and pathogenic mechanisms of group A Streptococcus.** *Clin Microbiol Rev.* 2014;27(2):264–301.
3. Thompson TZ, McMullen AR. **Group a streptococcus testing in pediatrics: The move to point-of-care molecular testing.** *J Clin Microbiol.* 2020;58(6).
4. Cohen JF, Bertille N, Cohen R, Chalumeau M. **Rapid antigen detection test for group A streptococcus in children with pharyngitis.** *Cochrane Database Syst Rev.* 2016;7(7).
5. Bessen DE, McShan WM, Nguyen SV, Shetty A, Agrawal S TH. **Molecular Epidemiology and Genomics of Group A Streptococcus.** *Infect Genet Evol.* 2015;33:393–418.

12. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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