

# V-Quick32

## Automatic nucleic acid extraction instrument



Integrated touchscreen.  
**User friendly** software



DNA/RNA **Universal Extraction Protocol**




**Fast extraction protocol:**  
Only 27 minutes



**Highly efficient:**  
Based on magnetic  
Beads technology



**Flexible Capacity:**  
From 1 to 32 samples



## VIASURE Universal Pathogens Extraction Kit for VIASURE V-Quick32, 96 prep

Ref. VS-EAU0196VQ



**Ready to use reagents**  
with preloaded plates.



**Compatible** with a wide  
range of sample types.



**Stable at Room  
Temperature.**



Purification of viral, bacterial,  
parasitic and fungal  
**DNA and RNA.**

